

# Anti-EBV Mixed Titer Performance Panel PME202

## A SERACARE PANEL PRODUCT

### INTENDED USE

Anti-EBV Mixed Titer Performance Panel PME202 is a panel of naturally occurring plasma specimens. PME202 is intended for use by diagnostic manufacturers and clinical laboratories to evaluate their Epstein-Barr Virus (EBV) test systems using well characterized specimens, and to provide comprehensive data for comparative analysis. For Research Use Only. Not for use in diagnostic procedures.

### PRODUCT DESCRIPTION

PME202 consists of a set of 21 aliquots from plasma units with reactivity ranging from negative to strongly positive for anti-EBV markers. Specimens are undiluted aliquots from plasma units. Units were aseptically filtered; no preservative was added.

Cat. No. PME202-1.2 1 vial per member  
21 members, 1.2 mL per vial

### STORAGE

Panel members should be stored frozen at -20°C or colder. SeraCare Life Sciences recommends that the panel members be divided into smaller aliquots to avoid multiple freeze-thaw cycles, if appropriate. If turbidity or particulate matter is observed, the samples should be centrifuged in accordance with each test kit manufacturer's instructions for sample preparation.

### INTERPRETATION OF RESULTS

The Data Sheet for anti-EBV Mixed Titer Performance Panel PME202-1.2 is available at [www.seracarepanels.com](http://www.seracarepanels.com). Select PME202 from the list of panel data sheets. The Data Sheet lists results for panel

members generated using commercially available screening and monitoring tests approved in the U.S. and/or the European Union. The tests that are listed were performed at SeraCare or at recognized reference laboratories (RL) by individuals who routinely use these procedures. Information regarding specific test methods is available on the Data Sheet.

### LIMITATIONS

PME202 is offered for research use only. Not for use in diagnostic procedures. Data are provided for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.

### PRECAUTIONS

Follow Universal Precautions<sup>1</sup>. The units that make up this panel were tested and found negative for HBsAg, anti-HIV 1/2, and anti-HCV. This does not ensure the absence of these or other human pathogens.

Do not pipette by mouth. Do not smoke, eat, or drink in areas where specimens are handled.

These materials should be disposed of in a manner that will inactivate pathogenic agents.

### REFERENCES

1. CDC recommendations for prevention of HIV transmission in health care settings. MMWR 36 (supp.2) 1987.

**For assistance, contact SeraCare Technical Support at 508.244.6400.**

The Data Sheet for this panel in PDF form can be found at [www.seracarepanels.com](http://www.seracarepanels.com)

The printed Data Sheet may be requested by email at [info@seracare.com](mailto:info@seracare.com), or by phone at 508.244.6400

### ASK ABOUT RELATED SERACARE PRODUCTS

- ACCURUN<sup>®</sup> independent quality controls
- SeraCon<sup>™</sup> and other processed plasma products
- Global Patient Sample (GPS) Program – access to a vast and evolving inventory of single test patient samples

