

# HIV-1 INCIDENCE/PREVALENCE PERFORMANCE PANEL PRB601

**BBI Diagnostics**  
*A Boston Biomedica Company*  
TEL: (508) 580-1900  
US Toll Free: (800) 676-1881  
FAX: (508) 580-2202

## INTENDED USE

The HIV-1 Incidence/Prevalence Performance Panel PRB601, is a panel of naturally occurring plasma specimens. This Panel is intended for use by diagnostic manufacturers and clinical laboratorians to evaluate assays that determine incidence (new infections) vs. prevalence (long-standing infections) using well-characterized specimens, and to provide comprehensive data for comparative analysis. For Research Use Only. Not for use in diagnostic procedures.

## PRODUCT DESCRIPTION

PRB601 consists of a set of 15 aliquots of plasma from different donors: seven members are characterized as incident and eight as prevalent, based on consensus results from nine tests using five different methods. Specimens are undiluted aliquots from plasma units collected from HIV positive deferred plasma donors in the United States whose dates of infection and seroconversion are unknown. Units were aseptically filtered and no preservative was added.

## REAGENTS

Cat. No. PRB601-1.0      1 vial per member  
   15 members, 1.0 ml per vial

## PRECAUTIONS

These materials have not been inactivated and should be considered biohazardous. Use the Centers for Disease Control and Prevention (CDC) recommended universal precautions for handling this product.<sup>1</sup> Panel members were tested and all were found positive by tests for anti-HIV 1.

The units that make up this panel were tested for HBsAg and HCV. All panel members were found to be non-reactive for HBsAg. One HCV reactive member is indicated on the data sheet.

Do not pipette by mouth.

Do not smoke eat or drink in areas where specimens are handled.

These materials should be disposed of in a manner

that will inactivate pathogenic agents.

## STORAGE

Panel members should be stored frozen at  $-20^{\circ}\text{C}$  or colder. BBI recommends that the panel members be divided into smaller aliquots for ease in repeated testing and to avoid multiple freeze-thaw cycles. If turbidity or particulate matter is observed, the samples should be centrifuged in accordance with each test kit manufacturer's instructions for sample preparation.

## INTERPRETATION OF RESULTS

The attached data sheet (page1) provides the results of assays for HIV-1 antigen, HIV RNA and HIV antibody generated using commercially available screening and confirmatory tests. Pages 2 and 3 provide results from incidence/prevalence tests. Testing was performed at BBI and at internationally recognized non-commercial referee laboratories (RL) by laboratorians who routinely use these test procedures

## LIMITATIONS

This Panel is offered for Research Use Only. Not for use in diagnostic procedures. The data are provided for informational purposes only. BBI Diagnostics does not claim that others can duplicate these results exactly.

## REFERENCES

1. CDC recommendations for prevention of HIV transmission in health care settings. MMWR 36 (sup.2) 1997.

For assistance, contact BBI Technical Support



375 West Street, West Bridgewater, MA 02379 USA  
E-mail: info@bbii.com, website: www.bbii.com

June 2004  
11122-01

**Boston Biomedica, Inc.**  
**HIV-1 INCIDENCE/PREVALENCE PERFORMANCE PANEL (PRB601)**  
**DATA SHEET**

**HIV DIAGNOSTIC TESTS**

-----US FDA Licensed HIV Tests-----

-----results expressed as s/co-----

Member I.D.#	Abbott	Bio-Rad	Coulter	Calypte	Roche
	HIV 1/2	HIV 1/2	HIV Ag	Anti-HIV-1 Western blot	HIV-1 RNA
	<u>BBi</u>	<u>BBi</u>	<u>BBi</u>	<u>BBi</u>	c/ml
					<u>BBi</u>
PRB601-01	18.5	>9.9	0.2	18,24,31,41,51,55,65,120,160	3 x 10 <sup>4</sup>
PRB601-02	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	4 x 10 <sup>3</sup>
PRB601-03	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	<400
PRB601-04	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	4 x 10 <sup>3</sup>
PRB601-05	15.8	>9.9	1.8	18,24,31,41,51,55,65,120,160	4 x 10 <sup>5</sup>
PRB601-06	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	2 x 10 <sup>4</sup>
PRB601-07	18.5	>9.9	1.1	18,24,31,41,51,55,65,120,160	2 x 10 <sup>5</sup>
PRB601-08	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	1 x 10 <sup>5</sup>
PRB601-09	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	1 x 10 <sup>4</sup>
PRB601-10	17.4	>9.9	0.7	18,24,31,41,51,55,65,120,160	4 x 10 <sup>4</sup>
PRB601-11	18.5	>9.9	0.2	18,24,31,41,51,65,120,160	5 x 10 <sup>3</sup>
PRB601-12	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	2 x 10 <sup>3</sup>
PRB601-13	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	3 x 10 <sup>3</sup>
PRB601-14	17.4	>9.9	1.3	18,24,31,41,51,55,65,120,160	3 x 10 <sup>5</sup>
PRB601-15	17.4	>9.9	0.2	18,24,31,41,51,55,65,120,160	4 x 10 <sup>4</sup>

  

<b>Run Date:</b>	6-Mar-04	4-Mar-04	12-Mar-04	18-Mar-04
<b>Kit Lot#:</b>	10905M100	318WBB-05	2063K194	A4316-394
<b>Exp.Date:</b>	15-Apr-04	13-Aug-04	16-Jun-04	11-Nov-04
<b>Product #:</b>	3A77-28	32588	6604535	98002

**CAUTION:** Panel members were tested and were found positive by tests for HIV. Panel member 8 was found reactive for HCV.

For research use only. Not for use in diagnostic procedures. Data are offered for informational purposes. BBI does not claim that others can duplicate these tests results exactly. Testing was performed at BBI or at internationally recognized referee laboratories (RL) by individuals who routinely use these procedures. EIA results are expressed as specimen absorbance to cutoff ratios (s/co). Ratios  $\geq 1.0$  are considered reactive. Specimens are undiluted aliquots from plasma units collected from 1996 through 1998.

## HIV INCIDENCE TESTS

		-----LESS SENSITIVE (LS)-----				-----AVIDITY INDEX (AI)-----		-----OTHER APPROACHES-----		
Member	Mfg.: Method: Units:	Abbott 3A11 SOD HIV-1	BioMerieux Vironostika SOD HIV-1	BioMerieux Vironostika SOD HIV-1/2	BioMerieux Vironostika+O SOD HIV-1/2	Bio-Rad Genetic Systems % rLAV	BioMerieux Vironostika % HIV-1/2	HIV RI Mean OD HIV 1	In-House BED CEIA Mean OD-n HIV 1	NRL µg/ml HIV-1
<u>I.D.#</u>	<u>CONSENSUS</u>	<u>RL11</u>	<u>RL7</u>	<u>RL36</u>	<u>RL 20</u>	<u>RL7</u>	<u>RL20</u>	<u>RL38</u>	<u>RL20</u>	<u>RL36</u>
PRB601-01	I	0.156	0.180	0.166	0.365	24	12	0.220	0.364	24
PRB601-02	I	0.051	0.100	0.093	0.195	12	18	0.160	0.132	36
PRB601-03	P	1.204	2.240	1.616	2.710	89	80	0.680	1.091	0
PRB601-04	P	1.163	1.910	2.148	4.459	94	80	0.970	1.412	0
PRB601-05	I	-0.060	0.010	0.028	0.187	10	8	0.160	0.332	67
PRB601-06	P	1.314	1.840	1.436	1.882	77	45	0.870	1.171	13
PRB601-07	I	0.209	0.220	0.186	0.294	30	20	0.260	0.239	28
PRB601-08	P	5.649	17.630	5.893	8.120	97	87	0.990	2.994	0
PRB601-09	I	0.349	0.340	0.261	0.526	46	18	0.340	0.436	35
PRB601-10	P	3.172	7.110	4.881	8.653	98	93	0.990	3.078	2
PRB601-11	P	2.376	4.790	2.793	5.015	86	77	0.900	1.213	11
PRB601-12	I	0.200	0.090	0.166	0.523	36	27	0.280	0.345	41
PRB601-13	P	5.100	11.080	5.512	6.785	98	90	0.990	3.034	0
PRB601-14	I	0.153	0.010	0.079	0.206	25	10	0.200	0.743	32
PRB601-15	P	1.021	1.920	1.566	2.544	92	72	0.940	1.683	0
<b>CUTOFF</b>		0.75	1.00	1.00	1.00	55	30	0.500	1.00	20

### DESCRIPTION OF METHODS

**Less Sensitive (LS)** tests involve "detuning" a commercially available HIV method to make it less sensitive by diluting samples, so that samples from early infection (positive on an unmodified test) will be below a specified cutoff, while samples from a long-standing infection will be above the cutoff. Results are expressed as Standard Optical Density (SOD) using a calibrator to normalize results.<sup>1,2</sup>

**Avidity Index (AI)** tests exploit the lower avidity of early vs. "mature" antibodies for the target HIV protein antigens. In AI methods, the sample is split, one aliquot is treated with a chaotropic agent in buffer and the other with buffer alone. Both aliquots are tested using a commercially available anti-HIV test method, and results are expressed as a percent of treated to untreated sample. In early infection samples, antibody-antigen complexes are disrupted by the chaotropic agent.<sup>3,4</sup>

**HIV RI (recent infection)** is an EIA that detects antibody responses to HIV epitopes empirically found to be expressed in early infection more intensely than in long-standing infection. In the HIV RI method, results are expressed as mean OD and compared to a calibrator.<sup>5</sup>

**BED-CEIA** is an IgG-capture EIA in which increased amounts of gp41-IDR specific antibodies are detected by BED-biotin peptide. Results are normalized using a calibrator OD and expressed as OD-n.<sup>6</sup>

**NRL HIV-1** is an assay that identifies individuals infected with HIV-1 within 120+/-16 days by detecting a specific antigen-antibody interaction that is transient in recent infection. Results are quantified using a standard curve, and values above 20 µg/ml are considered indicative of early infection. To identify HIV-1 infected individuals prior to day 30, this assay must be performed in parallel with a p24 antigen assay.<sup>7</sup>

## HIV INCIDENCE TESTS

Member I.D.#	-----LESS SENSITIVE (LS)-----				-----AVIDITY INDEX (AI)-----		-----OTHER APPROACHES-----			
	Mfg.: Method: Units:	Abbott 3A11 SOD HIV-1 RL11	BioMerieux Vironostika SOD HIV-1 RL7	BioMerieux Vironostika SOD HIV-1/2 RL36	BioMerieux Vironostika+O SOD HIV-1/2 RL 20	Bio-Rad Genetic Systems % rLAV RL7	BioMerieux Vironostika % HIV-1/2 RL20	HIV RI Mean OD HIV 1 RL38	In-House BED CEIA Mean OD-n HIV 1 RL20	NRL µg/ml HIV-1 RL36
<b>INCIDENT</b>										
PRB601-01		0.156	0.180	0.166	0.365	24	12	0.220	0.364	24
PRB601-02		0.051	0.100	0.093	0.195	12	18	0.160	0.132	36
PRB601-05		-0.060	0.010	0.028	0.187	10	8	0.160	0.332	67
PRB601-07		0.209	0.220	0.186	0.294	30	20	0.260	0.239	28
PRB601-09		0.349	0.340	0.261	0.526	46	18	0.340	0.436	35
PRB601-12		0.200	0.090	0.166	0.523	36	27	0.280	0.345	41
PRB601-14		0.153	0.010	0.079	0.206	25	10	0.200	0.743	32
<b>PREVALENT</b>										
PRB601-03		1.204	2.240	1.616	2.710	89	80	0.680	1.091	0
PRB601-04		1.163	1.910	2.148	4.459	94	80	0.970	1.412	0
PRB601-06		1.314	1.840	1.436	1.882	77	45	0.870	1.171	13
PRB601-08		5.649	17.630	5.893	8.120	97	87	0.990	2.994	0
PRB601-10		3.172	7.110	4.881	8.653	98	93	0.990	3.078	2
PRB601-11		2.376	4.790	2.793	5.015	86	77	0.900	1.213	11
PRB601-13		5.100	11.080	5.512	6.785	98	90	0.990	3.034	0
PRB601-15		1.021	1.920	1.566	2.544	92	72	0.940	1.683	0
<b>CUTOFF</b>		0.75	1.00	1.00	1.00	55	30	0.500	1.00	20

## References:

- Janssen RS, Satten GA, Stramer SL, Rawal BD, O'Brien TR, Weiblen BJ, Hecht FM, Jack N, Cleghorn FR, Kahn JO, Chesney MA, Busch MP. New testing strategy to detect early HIV-1 infection for use in incidence estimates and for clinical and prevention purposes. *JAMA* 1998 Jul 1;280(1):42-8: Erratum: *JAMA* 1999 May 26;281(20):1893.
- Kothe D, Byers RH, Caudill SP, Satten GA, Janssen RS, Hannon WH, Mei JV. Performance Characteristics of a New Less Sensitive HIV-1 Enzyme Immunoassay for Use in Estimating HIV Seroincidence. *JAIDS* 2003 Aug 15;33(5):625-34.
- Suligoi B, Galli C, Massi M, Di Sora F, Sciandra M, Pezzotti P, Recchia O, Montella F, Sinicco A, Rezza G. Precision and Accuracy of a Procedure for Detecting Recent Human Immunodeficiency Virus Infections by Calculating the Antibody Avidity Index by an Automated Immunoassay-Based Method. *JCM* 2002 Nov;40(11):4015-20.
- Jenner J, Grazioplene M, Kazianis A, Phinney K, and Werner B. Modification of a Commercial HIV-1 EIA for Identification of Recent HIV-1 Infections: Use of Differential Antibody Avidity. Poster presentation at 10th Conference on Retroviruses and Opportunistic Infections, Boston, MA, February 10-14, 2003.
- Barin F, Meyer L, Lancar R, Deveau C, Gharib M, Laporte A, Desenclos JC, Costagliola D. A New Immunoassay for the Identification of Recent HIV-1 Infections: Development and Validation. Poster presentation LB21 at The 2nd IAS Conference on HIV Pathogenesis and Treatment, Paris, July 13-16, 2003.
- Hu DJ, Vanichseni S, Mock PA, Young NL, Dobbs T, Byers R, Choopanya K, Griensven F, Kitayaporn D, McDougal JS, Tappero JW, Mastro TD, Parekh BS. HIV-1 incidence estimates by detection of recent infection from a cross-sectional sampling of injection drug users in Bangkok: Use of the IgG Capture BED enzyme immunoassay. *AIDS Res and Hum Retro* 2003;19:727-730.
- Wilson KM, Croom HA, Richards K, Doughty L, Cunningham PH, Kemp BE, Branson BM, Johnson EIM, Dax EM. Incidence Immunoassay for Distinguishing Recent from Established HIV-1 Infection in Untreated Populations. (submitted for publication)