

NOTEBOOK-READY PROCEDURE

Sandwich ELISA Procedure

Experiment/Project #: _____

Name(s): _____

Date(s): _____

Material	Company/Catalog #/Lot #
Plate:	
Coating Antibody:	
Antigen/Sample:	
Secondary Antibody:	
Tertiary Antibody:	
Wash Solution:	
Blocker:	
Substrate:	
Coating Solution:	
Antibody Diluent:	
Stop Solution:	
Equipment	Company/Catalog #/Lot #
Plate Washer:	
Plate Reader:	

1. **DETECTION SYSTEM:** Chemiluminescent Colorimetric Fluorescent

2. **PLATE TYPE (96 WELL)**

Clear Black White
 Flat Bottom U Bottom C Bottom
 High Binding Medium Binding Other: _____

3. **ANTIBODY COATING**

1 hour at room temperature (RT) or overnight 4°C Temperature: _____
 Coating Antibody: _____ Concentration: _____
 Coating Solution: _____ Volume/Well: _____
 Antibody: _____ Polyclonal Monoclonal
 Antibody Purity: Unpurified Serum Protein A/G Antigen Affinity Ascites
 Actual time: Start - _____ AM/PM End - _____ AM/PM Total - _____

4. **BLOCK PLATE**

1 hour at RT or overnight 4°C Blocker: _____ Temperature: _____
 Volume/Well: _____
 Actual time: Start - _____ AM/PM End - _____ AM/PM Total - _____

5. **ANTIGEN/SAMPLE INCUBATION**

30 – 60 minutes at RT or overnight 4°C Temperature: _____
 Antigen: _____ Diluent: _____
 Antigen Concentration: _____ Volume/Well: _____
 Actual time: Start - _____ AM/PM End - _____ AM/PM Total - _____

6. WASH PLATE

Wash Solution: _____ Volume/Well: _____
Actual Time: Wash 1 - _____ Wash 2 - _____ Wash 3 - _____

7. PRIMARY ANTIBODY INCUBATION

30 – 60 minutes at RT or overnight 4°C Temperature: _____
Antibody: _____ Polyclonal Monoclonal
Antibody Diluent: _____
Conjugate: N/A HRP AP Other Tag: _____
Antibody Purity: Unpurified Serum Protein A/G Antigen Affinity Ascites
Antibody Concentration: _____ Dilution: _____
Actual Time: Start - _____ AM/PM End - _____ AM/PM Total - _____

8. WASH PLATE

N/A
Wash Solution: _____
Actual Time: Wash 1 - _____ Wash 2 - _____ Wash 3 - _____

9. SECONDARY ANTIBODY INCUBATION

30 – 60 minutes at RT or overnight 4°C Temperature: _____
Antibody: _____ Polyclonal Monoclonal
Antibody Diluent: _____
Conjugate: HRP AP Other Tag: _____
Antibody Purity: Unpurified Serum Protein A/G Antigen Affinity Ascites
Antibody Concentration: _____ Dilution: _____
Actual Time: Start - _____ AM/PM End - _____ AM/PM Total - _____

10. WASH PLATE

N/A
Wash Solution: _____
Actual Time: Wash 1 - _____ Wash 2 - _____ Wash 3 - _____

11. SUBSTRATE INCUBATION

Substrate Name: _____ Colorimetric Chemiluminescent Fluorescent
Recommended development time: _____ minutes Actual Time: _____ minutes
Stop Solution: _____ N/A

12. DATA ANALYSIS

Time: _____ minutes Slit Width: _____ Time/Well: _____
Computer File Name: _____ Location: _____
Wavelength: _____

13. RESULTS/CONCLUSIONS




www.seracare.com
800.676.1881
508.244.6400



Power Your
Immunoassays

www.kpl.com
800.638.3167
301.948.7755

14. PLATE TEMPLATE



www.seracare.com
800.676.1881
508.244.6400



Power Your
Immunoassays

www.kpl.com
800.638.3167
301.948.7755