

BOSTON BIOMEDICA, INC.
ANTI-HIV 1 SEROCONVERSION PANEL N (PRB914)
DATA SHEET

-----Test Kits Licensed by the U.S. FDA-----

Panel	Bleed	Days Since	Abbott HIV BBI	Abbott HIV 1/2 BBI	CBC HIV BBI	CPI HIV BBI	Gen. Sys. HIV BBI	Gen. Sys. HIV 1/2 BBI	Org. Tek. HIV BBI	Syva HIV BBI	Abbott HIV-Ag BBI
<u>I.D. #</u>	<u>Dates</u>	<u>1st Bleed</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>
PRB914-01	01/12/90	0	1.4	8.4	1.6	0.7	0.6	0.8	0.7	1.9	0.4
PRB914-02	01/16/90	4	2.0	7.8	1.8	1.1	0.5	1.1	0.9	1.9	0.5
PRB914-03	01/19/90	7	2.1	8.0	2.0	1.1	0.7	1.4	0.9	2.1	0.5
PRB914-04	02/06/90	25	5.3	12.8	2.8	3.2	2.1	2.8	2.2	2.9	0.4
PRB914-05	02/12/90	31	4.7	11.4	2.2	3.2	1.8	3.0	2.1	2.0	0.4
Run Date			4/13/92	2/27/92	1/10/92	1/13/92	1/17/92	11/5/91	12/28/91	1/10/92	5/10/91
Kit Lot #			61399M201	60147M101	A6857	J24713	258GX1	171GR1	120477	A6853	51907M201
Exp. Date			4/15/92	5/19/92	8/10/92	7/15/92	3/31/92	3/31/92	7/92	6/12/92	7/8/91

-----Test Kits Licensed in Europe-----

-----FDA-Licensed Confirmatory Tests-----

Panel	Behring HIV 1/2 RL	IAF Biochem HIV 1/2 BBI	Innotest HIV 1/2 RL	Wellcozyme HIV 1/2 RL	BioRad Western Blot BBI	Dupont Western Blots BBI
<u>I.D. #</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>s/co</u>	<u>Band Patterns</u> <u>Result</u>	<u>Band Patterns</u> <u>Result</u>
PRB914-01	3.9	4.4	1.7	5.1	f24 IND	f24 IND
PRB914-02	4.2	4.0	2.2	6.3	24, vf55 IND	24, vf160 POS
PRB914-03	4.0	3.8	2.1	7.2	24, f55 IND	24, vf160 POS
PRB914-04	7.4	10.7	3.9	9.2	18, 24, 55, vf160 POS	18, 24, f120, 160 POS
PRB914-05	7.5	10.4	3.8	9.5	18, 24, 55, vf160 POS	18, 24, 120, 160 POS
Run Date		Not	5/2/91	Not	Not	5/31/91
Kit Lot #		Available	1055-2	Available	Available	B1220-161
Exp. Date			2/19/92			8/8/91

Panel N has tested negative for HBsAg, anti-HCV and anti-HTLV.

Western Blots are interpreted using ASTPHLD/CDC criteria (MMWR, vol. 38, S-7, 1989).

vf = very faint, f = faint

For investigational use only. Data are offered for informational purposes. BBI does not claim that others can duplicate these test results exactly. EIA and Western Blot results were generated using commercially available FDA approved tests, performed at BBI and at an internationally recognized non-commercial Referee Laboratory (RL) by individuals who routinely use these test procedures. All numeric results are means of duplicates, expressed as specimen absorbance to cutoff ratios. Ratios ≥ 1.0 are considered reactive. Specimens are undiluted aliquots from plasma units collected from a single donor in 1990. No preservatives were added.