

# AccuVert™ HCV Seroconversion Panel

PHV925 (0610-0224) / Batch #10375695

## OVERVIEW


AccuVert™ HCV Seroconversion Panel PHV925 (0610-0224) / Batch #10375695 is a 5-member panel of undiluted, naturally occurring plasma samples (1 vial per member, 1.2 mL) collected over 27 days from a 24-year old in 2008 during a period of HCV seroconversion. This donor has been characterized as genotype 1a.

Test results from commercially-available HCV RNA, antibody and confirmatory assays are included for reference. This panel converts from negative to positive for anti-HCV.

For Research Use Only. Not for use in diagnostic procedures. Data are offered for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: Potentially infectious materials. Follow Universal Precautions. The units that make up this panel were tested and found negative for Anti-HIV-1/2 and HBsAg. This does not ensure the absence of these or other human pathogens.

Evolution of HCV Markers in Early Infection



This graph demonstrates the evolution of early HCV infection amongst panel members from the Abbott ARCHITECT (Anti-HCV) and Abbott m2000 (HCV RNA) test methods.

# AccuVert™ HCV Seroconversion Panel

PHV925 (0610-0224) / Batch #10375695

## Panel Member Information

| Panel Member | SeraCare Batch # | SeraCare Donor ID # | Bleed Date  | Days Since 1 <sup>st</sup> Bleed |
|--------------|------------------|---------------------|-------------|----------------------------------|
| 01           | 9239419          | BD107566            | 03-Dec-2008 | 0                                |
| 02           | 9239420          | BD107566            | 05-Dec-2008 | 2                                |
| 03           | 9239421          | BD107566            | 11-Dec-2008 | 8                                |
| 04           | 9239423          | BD107566            | 13-Dec-2008 | 10                               |
| 05           | 9239424          | BD107566            | 30-Dec-2008 | 27                               |

## HCV RNA and HCV Confirmatory

| Panel Member          | Abbott<br>m2000<br>Realtime HCV<br>(IU/mL) <sup>1,2</sup> | Roche COBAS®<br>AmpliPrep/COBAS®<br>TaqMan®<br>HCV Quantitative Test,<br>v2.0<br>(IU/mL) <sup>1,2</sup> | Fujirebio<br>INNO-LIA®<br>HCV Score<br>Band Pattern | Fujirebio<br>INNO-LIA®<br>HCV Score<br>Interpretation |
|-----------------------|---|---|---|---|
|                       |   | 01  | <b>8.6E+05</b>                                      | <b>1.0E+06</b>  |
| 02                    | <b>1.5E+06</b>  | <b>3.4E+06</b>  | No Bands  | NEG   |
| 03                    | <b>3.1E+06</b>  | <b>3.7E+06</b>  | No Bands  | NEG   |
| 04                    | <b>3.0E+06</b>  | <b>3.8E+06</b>  | No Bands  | NEG   |
| 05                    | <b>7.6E+05</b>  | <b>8.6E+05</b>  | <b>C1, C2, NS3, NS4, NS5</b>                        | <b>POS</b>  |
| Test Date             | 01-Nov-2018   | 01-Nov-2018   | 15-Nov-2018   |   |
| Test Site             | RL  | RL  | RL  |   |
| Kit Part Code         | NA  | NA  | NA  |   |
| Kit Lot No.           | NA  | NA  | 404975  |   |
| Kit Exp. Date         | NA  | NA  | 28-Feb-2019   |   |
| Kit Regulatory Status | IVD   | IVD/CE  | IVD/CE  |   |

<sup>1</sup>Results are reported as international units per mL (IU/mL); positive/reactive results are noted in bold red.

<sup>2</sup>Results are reported as the mean result of duplicate testing.

NA = Not Available

RL = Reference Lab

IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE Marking

AccuVert™ HCV  
Seroconversion Panel  
PHV925 (0610-0224) / Batch #10375695

## Anti-HCV

| Panel Member          | Abbott<br>ARCHITECT<br>Anti-HCV<br>(s/co) <sup>1,3</sup> | Ortho®<br>HCV<br>Version 3.0<br>ELISA Test System<br>(s/co) <sup>1,3</sup> | Roche<br>Elecsys®<br>Anti-HCV<br>(s/co) <sup>1,3</sup> | Siemens<br>ADVIA Centaur®<br>Anti-HCV<br>(INDEX) <sup>2,3</sup> |
|-----------------------|--|--|--|---|
| 01                    | 0.1  | 0.0  | 0.1  | 0.1   |
| 02                    | 0.1  | 0.0  | 0.1  | 0.1   |
| 03                    | 0.2  | 0.0  | <b>2.3</b>   | 0.1   |
| 04                    | 0.3  | 0.0  | <b>9.4</b>   | 0.3   |
| 05                    | <b>9.0</b>   | <b>&gt;6.5<sup>4</sup></b>   | <b>56.8</b>  | <b>&gt;11.0<sup>4</sup></b>                                     |
| Test Date             | 29-Oct-2018  | 26-Oct-2018  | 08-Nov-2018  | 31-Oct-2018   |
| Test Site             | SC   | SC   | RL   | RL  |
| Kit Part Code         | 1L79   | 930740   | NA   | NA  |
| Kit Lot No.           | 83370LI00  | TXE670   | NA   | NA  |
| Kit Exp. Date         | 15-Nov-2018  | 10-May-2019  | NA   | NA  |
| Kit Regulatory Status | IVD/CE   | IVD  | IVD  | IVD/CE  |

<sup>1</sup>Results are reported as a signal to cutoff ratio (s/co); positive/reactive results are noted in bold red.

<sup>2</sup>Results are reported as an Index value (INDEX); positive/reactive results are noted in bold red.

<sup>3</sup>Results are reported as the mean result of duplicate testing.

<sup>4</sup>Off-scale

NA = Not Available

SC = SeraCare; RL = Reference Lab

IVD = In Vitro Diagnostic; CE = Conformité Européenne or CE Marking

The package insert for this panel can be found at [www.seracare.com](http://www.seracare.com).

A printed copy of the package insert or data sheet may be requested by email at [info@seracare.com](mailto:info@seracare.com) or by phone at 508.244.6400.