

AccuSet™ HBsAg Performance Panel

0805-0391 / Batch #10494991

OVERVIEW


AccuSet™ HBsAg Performance Panel (0805-0391) is a 25-member validation panel of undiluted, naturally occurring plasma samples (1 vial per member, 1.5 mL per vial). Panel members represent bleeds from multiple individuals positive for HBsAg. Each sample represents a single collection event. No preservatives were added.

Test results from commercially-available HBV assays are included for characterization of the panel members. This panel of human plasma samples demonstrates a range of reactivity for HBsAg test methods. Additional tests performed on HBV DNA, HBeAg, Anti-HBc IgM, Anti-HBc and Anti-HBs are also provided. One sample is included as a non-reactive sample and is negative for all HBV test methods performed.

For Research Use Only. Not for use in diagnostic procedures. Data are provided for informational purposes. SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: Potentially infectious materials. Follow Universal Precautions. The units that make up this panel were tested and found negative for anti-HIV-1/2 and anti-HCV. This does not ensure the absence of these or other human pathogens.

AccuSet™ HBsAg Performance Panel



This graph demonstrates reactivity amongst panel members from the Roche Elecsys 2010 HBsAg test method.

AccuSet™ HBsAg
Performance Panel
0805-0391 / Batch #10494991

Panel Member Information

Panel Member	SeraCare Batch #	SeraCare Donor ID #	Collection Date
01	10266975	BD328487	14-Sep-2015
02	10267001	BD328495	18-Jan-2016
03	9107572	BD110467	26-Nov-1999
04	9107562	BD110465	NA
05	BM204966	BD102166	12-Jan-2006
06	9259534	BD111350	26-Nov-2010
07	10267027	BD328507	26-Dec-2015
08	9120448	BD110515	NA
09	BM206376	BD111334	30-Jun-2006
10	9121493	BD110525	NA
11	BM216724	NA	NA
12	9120182	BD110510	NA
13	9121507	BD110529	NA
14	9253281	BD110945	18-Feb-2011
15	9233766	BD107191	16-Oct-2001
16	10266974	BD328486	20-Mar-2016
17	10267015	BD328502	21-Nov-2015
18	BM217727	BD105187	18-Aug-2007
19	9118028	BD110506	NA
20	9124953	NA	NA
21	BM205270	BD102507	26-Apr-2006
22	9253283	BD110947	14-Jan-2011
23	9242009	BD110034	29-Jun-2009
24	10266961	BD328480	28-Oct-2015
25	10310183	BD349223	31-May-2016

NA = Not Available

AccuSet™ HBsAg Performance Panel

0805-0391 / Batch #10494991

HBV DNA and HBsAg

Panel Member	Roche COBAS® AmpliPrep/COBAS® TaqMan® HBV Test, v2.0 (IU/mL) ^{1,4}	Roche Elecys® 2010 HBsAg (s/co) ^{2,4}	Siemens ADVIA Centaur HBsAg (Index) ^{3,4}	Abbott ARCHITECT HBsAg Qualitative (s/co) ^{2,4}
01	9.78E+03	3719.0	> 1000⁵	6435.8
02	3.01E+05	3215.0	> 1000⁵	4608.9
03	4.74E+04	47.7	1.32E+02	97.5
04	5.23E+04	313.5	> 1000⁵	714.5
05	NA	2.0	3.92E+00	3.4
06	2.51E+02	5.2	1.53E+01	6.6
07	TND	126.0	5.20E+02	261.6
08	9.90E+03	34.0	9.91E+01	61.6
09	2.66E+04	4188.0	> 1000⁵	5529.6
10	2.61E+04	23.0	3.94E+01	44.5
11	8.60E+03	5.2	1.19E+01	8.3
12	3.71E+04	236.0	6.91E+02	435.4
13	2.84E+04	58.7	1.92E+02	112.0
14	1.01E+05	4238.0	> 1000⁵	5416.6
15	1.43E+04	4780.0	> 1000⁵	4417.3
16	1.01E+05	4072.0	> 1000⁵	6111.4
17	2.95E+05	4059.0	> 1000⁵	5822.9
18	2.76E+05	44.6	1.61E+02	108.0
19	3.07E+03	7.1	1.74E+01	13.3
20	1.95E+04	22.6	6.76E+01	41.1
21	4.57E+04	78.5	2.26E+02	137.4
22	TND	1.2	1.92E+00	1.9
23	5.67E+02	1.5	2.70E+00	1.9
24	TND	5.4	1.20E+01	8.3
25	TND	0.4	1.40E-01	0.2
Test Date	01-Jun to 03-Jun-2020	22-Jun-2020	01-Jun-2020	01-Jun-2020
Test Site	RL	RL	RL	SC
Kit Part Code	NA	NA	NA	4P53
Kit Lot No.	F22485	NA	109215	11379FN00
Kit Exp. Date	31-Jan-2021	NA	09-Jan-2021	11-Aug-2020

¹Results are reported as International Units per milliliter (IU/mL); positive/reactive results are noted in bold red.

²Results are reported as signal-to-cutoff ratio (s/co); positive/reactive results are noted in bold red.

³Results are reported as Index; positive/reactive results are noted in bold red.

⁴Results are reported as the mean result of duplicate testing.

⁵Results are off scale.

TND = Target Not Detected, NA = Not Available

RL = Reference Lab; SC = SeraCare

AccuSet™ HBsAg Performance Panel

0805-0391 / Batch #10494991

HBsAg

Panel Member	Ortho VITROS® HBsAg (s/co) ^{1,2}	Bio-Rad Genetic Systems™ HBsAg EIA 3.0 (s/co) ^{1,2}
01	5575.0	>42.1³
02	2960.0	>42.1³
03	44.3	29.6
04	317.0	>42.1³
05	1.3	3.0
06	6.2	4.8
07	73.0	22.7
08	29.5	27.2
09	4025.0	>42.1³
10	12.2	16.3
11	3.0	5.3
12	217.5	>42.1³
13	42.1	27.6
14	4080.0	>42.1³
15	3035.0	>42.1³
16	5160.0	>42.1³
17	4800.0	>42.1³
18	85.0	39.5
19	5.4	7.6
20	19.4	24.2
21	71.8	39.3
22	0.2	6.3
23	0.7	12.7
24	4.2	12.2
25	0.1	0.3
Test Date	02-Jun-2020	29-May to 08-Jun-2020
Test Site	RL	SC
Kit Part Code	NA	32591
Kit Lot No.	NA	170PCC-05
Kit Exp. Date	NA	31-Aug-2020

¹Results are reported as signal-to-cutoff ratio (s/co); positive/reactive results are noted in bold red.

²Results are reported as the mean result of duplicate testing.

³Results are off scale.

POS = Positive; NEG = Negative

RL = Reference Lab; SC = SeraCare; NA = Not Available

AccuSet™ HBsAg Performance Panel

0805-0391 / Batch #10494991

HBeAg and Anti-HBe

Panel Member	Abbott ARCHITECT HBeAg (s/co) ^{1,3}	DiaSorin ETI-EBK Plus HBeAg EIA (s/co) ^{1,3}	DiaSorin ETI-AB-EBK Plus Anti-HBe EIA (co/s) ^{2,3}
01	0.25	0.02	225.00
02	132.68	65.57 ⁴	0.17
03	0.86	1.22	0.46
04	6.65	11.56	0.42
05	0.27	0.05	0.64
06	0.30	0.08	0.69
07	0.23	0.02	123.75
08	0.38	0.33	0.60
09	1.38	1.82	1.19
10	0.75	1.09	0.63
11	0.32	0.16	0.52
12	1.84	3.21	0.52
13	0.77	1.66	0.55
14	0.85	3.26	1.25
15	0.96	1.67	1.17
16	1.66	3.07	1.37
17	1.43	1.00 ⁵	2.17
18	15.79	28.80	0.40
19	0.32	0.12	0.48
20	0.69	1.33	0.49
21	2.33	3.61	0.59
22	0.25	0.02	0.52
23	0.25	0.02	0.51
24	0.23	0.02	6.93
25	0.26	0.02	0.43
Test Date	07-Jun-2020	29-May-2020 to 08-Jun-2020	01-Jun-2020
Test Site	RL	SC	SC
Kit Part Code	NA	P001930	P001929
Kit Lot No.	12002BE00	0600700A	0490790AA
Kit Exp. Date	10-Nov-2020	30-Sep-2020	16-Sep-2020

¹Results are reported as signal-to-cutoff ratio (s/co); positive/reactive results are noted in bold red.

²Results are reported as cutoff-to-signal ratio (co/s); positive/reactive results are noted in bold red.

³Results are reported as the mean result of duplicate testing.

⁴Result is off scale.

⁵Result is the average of four replicates (range 0.7-1.5 s/co).

RL = Reference Lab; SC = SeraCare; NA = Not Available

AccuSet™ HBsAg Performance Panel

0805-0391 / Batch #10494991

Anti-HBc IgM, Anti-HBc and Anti-HBs

Panel Member	DiaSorin	DiaSorin	DiaSorin
	ETI-CORE-IgMK PLUS Anti-HBc IgM EIA (s/co) ^{1,3}	ETI-AB-COREK PLUS Anti-HBc EIA (co/s) ^{2,3}	ETI-AB-AUK Plus Anti-HBs EIA (s/co) ^{1,3}
01	0.20	321.00	0.02
02	0.30	233.45	0.01
03	0.16	0.33	0.02
04	0.19	0.33	0.04
05	0.16	0.34	0.02
06	0.11	0.48	0.01
07	0.16	428.00 ⁴	0.02
08	0.16	0.45	0.02
09	0.36	428.00 ⁴	0.02
10	0.15	0.40	0.01
11	0.17	0.35	0.02
12	0.20	0.33	0.04
13	0.14	0.39	0.04
14	1.17	428.00 ⁴	0.01
15	1.11	428.00 ⁴	0.03
16	0.51	428.00 ⁴	0.01
17	2.52	428.00 ⁴	0.02
18	0.16	0.41	0.02
19	0.15	0.39	0.02
20	0.16	0.37	0.02
21	1.13	20.00	0.05
22	0.15	0.46	0.02
23	0.14	0.72	0.01
24	0.11	428.00 ⁴	0.01
25	0.19	0.48	0.06
Test Date	09-Jun-2020	04-Jun-2020	08-Jun-2020
Test Site	SC	SC	SC
Kit Part Code	P001928	P001927	P001931
Kit Lot No.	0270510A	9910710A	0650500A
Kit Exp. Date	30-Jun-2020	30-Sep-2020	30-Sep-2020

¹Results are reported as signal-to-cutoff ratio (s/co); positive/reactive results are noted in bold red.²Results are reported as cutoff-to-signal ratio (c/s); positive/reactive results are noted in bold red.³Results are reported as the mean result of duplicate testing.⁴Result is off scale

SC = SeraCare

The package insert for this panel can be found at www.seracare.com.

A printed copy of the package insert or data sheet may be requested by email at info@seracare.com or by phone at 508.244.6400.