# **sera** care

### **KPL DAPI**

(4', 6-diamidine-2'-phenylindole, dihydrochloride)

#### **DESCRIPTION**

KPL DAPI (4',6-diamidine-2'-phenylindole, dihydrochloride) is a fluorescent counterstain which stains the nucleus blue in color. The bright blue fluorescence provides excellent contrast to red fluorochromes such as rhodamine, Cy3, Texas Red, or R-Phycoerythrin.

#### FORM/STORAGE

KPL DAPI is supplied as a yellow-orange powder. Store powder at 2-8°C. In powdered form, KPL DAPI is stable for a minimum of 1 year from date of receipt.

#### RECONSTITUTION

Reconstitute KPL DAPI powder with 1.0 mL of reagent quality water to a final concentration of 1.0 mg/mL. The reconstituted stock solution should be stored at -20°C and is stable for approximately 6 months.

#### **APPLICATIONS**

KPL DAPI can be used for fluorescent *in situ* hybridization, immunohistochemistry, and flow cytometry.

#### SUGGESTED PROTOCOL

The dilutions and staining times stated here are offered as suggestions. Due to individual preference for stain intensity, these parameters should be optimized by the user.

#### Fluorescent In situ

#### Hybridization/Immunohistochemistry:

- Detect sample using an appropriate fluorescent conjugate.
- 2. Dilute KPL DAPI stock solution 1:1000 with reagent quality water, (i.e. 1.0  $\mu$ L DAPI + 999  $\mu$ L reagent quality water).
- Apply diluted KPL DAPI generously to the slide and incubate at room temperature for 3-5 minutes in the dark (cover slide with aluminum foil).
- 4. Rinse the slide for 5-10 seconds in 1X PBS.
- Mount with KPL Fluorescent Mounting Media or other appropriate mounting media.

## SOLUTION PREPARATION PBS:

KPL Propidium Iodide

To 500 mL reagent quality water, add 8.0 g NaCl, 0.2 g KCl, 1.44 g Na<sub>2</sub>HPO<sub>4</sub>, and 0.24 g KH<sub>2</sub>PO<sub>4</sub>. Adjust the pH to 7.4 with HCl. QS to 1 Liter with reagent quality water.

#### PRODUCT SAFETY AND HANDLING

This product is considered non-hazardous as defined by the Hazard Communication Standard (29 CFR 1910.1200).

Avoid contact with skin and eyes. In case of contact or spillage, flush thoroughly with copious amounts of water.

| RELATED PRODUCTS                             | CAT. NO.                 |
|--|--------------------------|
| KPL Fluorescent<br>Mounting Media            | 5570-0005 (71-00-16)     |
| KPL Cy3-Labeled<br>Streptavidin              | 5270-0022 (072-01-30-00) |
| KPL R-Phycoerythrin-<br>Labeled Streptavidin | 5270-0021 (0718-30-00)   |

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.

5930-0008 (71-04-01)

508.244.6400 • 800.676.1881 Toll Free • 508.634.3334 Fax

www.seracare.com