

# AccuPlex™ SARS-CoV-2 Quality Solutions

## ASSAY VERIFICATION AND ONGOING PERFORMANCE MONITORING FOR ROUTINE PATIENT TESTING

AccuPlex™ SARS-CoV-2 Verification Panel, coupled with the AccuPlex™ SARS-CoV-2 Molecular Controls Kit, is a complete Quality Solution designed to challenge the entire molecular test procedure ensuring clinical laboratories can have confidence in their SARS-CoV-2 assay results.

AccuPlex™ SARS-CoV-2 Verification Panel is optimized for assay verification at installation by documenting test performance along the assay's range enabling laboratories to establish lower limits of detection, perform assay comparisons, and evaluate staff proficiency. The product contains positive materials including the full SARS-CoV-2 viral genome, and negative materials targeting the human RNase P gene.

Product Description	Material Number	Sequence Coverage	Pack Size		Concentration
AccuPlex™ SARS-CoV-2 Verification Panel	0505-0168	Full Genome	Positive 1	1 x 3 mL	100,000 copies/mL
			Positive 2	1 x 3 mL	10,000 copies/mL
			Positive 3	1 x 3 mL	1,000 copies/mL
			Negative	1 x 3 mL	5,000 copies/mL (RNase P)

Not for In Vitro Diagnostic Use. Research Use Only.

AccuPlex<sup>™</sup> SARS-CoV-2 Molecular Controls Kit is designed to measure day-to-day performance of the assay. The product contains positive materials including the full SARS-CoV-2 viral genome, and negative materials targeting the human RNase P gene. The controls are formulated for use with in vitro diagnostic test methods that detect SARS-CoV-2 virus, the causative agent of COVID-19 disease.

Product Description	Material Number	Sequence Coverage	Pac	k Size	Concentration
AccuPlex™ SARS-CoV-2 Molecular Controls Kit	0505-0159	Full Genome			5,000 copies/mL 5,000 copies/mL (RNase P)

For In Vitro Diagnostic Use. CE-IVD Marked.

## HIGHLIGHTS

OFFERING FULL VIRAL
GENOME COVERAGE COMPATIBLE WITH ANY
SARS-COV-2 TARGET OF
INTEREST

REAL MAMMALIAN
VIRUS ENABLES FULL
PROCESS QC FOR RTPCR BASED VIRAL
DISEASE ASSAYS

NON-INFECTIOUS
AND REPLICATION
DEFICIENT; ENSURES
SAFETY DURING
ASSAY DEVELOPMENT
AND END-USER QC
APPLICATIONS

## **SARS-COV-2 ACCUPLEX SOLUTION:**

- · Non-infectious and replication deficient, enables safe handling in contrast to viral samples
- Fully-extractable with a real viral protein coat; superior to "naked" transcribed RNA
- 2 year stability at 2 8°C

To learn more about AccuPlex SARS-CoV-2 Quality Solutions visit: www.seracare.com/SARSCoV2 or contact customer service at 800.676.1881

your diagnostic assays. Compatible with multiplexed RT-PCR and NGS-based assays, AccuPlex custom recombinant materials are constructed with a replication-deficient mammalian virus, producing a safe, non-infectious material (Figure 1). With a protein coat and lipid bilayer, these mammalian virus-based reference materials resemble the complexity of virus targets found in true patient samples.

LGC SeraCare's AccuPlex recombinant material serves as a true full process molecular control for

# PATHOGENIC RNA VIRUS RNA VIRUS RNA Sequences 1. Target Virus Recombinant RNA Recombinant RNA Genome Target Virus RNA Sequences

FIGURE 1: 1) RNA sequence from the pathogenic virus of interest is chosen. 2) DNA synthesis and cloning occur to produce the recombinant RNA. 3) Recombinant RNA and helper RNA are co-transfected into the mammalian cells, allowing the encapsulation of recombined RNA. 4) Exocytosis of the mature enveloped non-infectious and replication deficient RNA virus containing the assay target RNA sequence of interest

## MOST 'PATIENT SAMPLE-LIKE' MATERIAL EVALUATES ENTIRE WORKFLOW

Unlike RNA-spiked buffer, or technologies that package viral RNA into a bacteriophage, the AccuPlex recombinant closely resembles the wild-type mammalian pathogenic virus. This enables the release of the viral genome at a similar rate to the wild-type virus during the nucleic acid sample preparation process. AccuPlex recombinant material mimics a real patient sample in your workflow, serving as a full-process control for your assay.

## ACCUPLEX SOLVES ASSAY DEVELOPMENT CHALLENGES

If you're developing diagnostics for emerging viral diseases and have the challenge of including safe, non-infectious controls in your test kit, partner with LGC SeraCare's talented R&D team to produce your custom AccuPlex recombinant virus material (DNA or RNA-based). Utilizing your sequences of interest and product specifications, we will develop a custom solution which meets your unique requirements.

## ORDERING INFORMATION

To place an order, or learn more about our SARS-CoV-2 Quality Solutions, please contact us at +1.508.244.6400 and 800.676.1881 or email <a href="mailto:customerservice@seracare.com">customerservice@seracare.com</a>



AccuPlex<sup>™</sup> is a trademark of SeraCare Life Sciences, Inc. © 2021 SeraCare Life Sciences, Inc. All rights reserved. MKT-00520-07

## ABOUT LGC SERACARE

TRUSTED SUPPLIER
TO THE DIAGNOSTIC
TESTING INDUSTRY
FOR OVER 30 YEARS

HIGH-QUALITY

CONTROL PRODUCTS,

RAW BIOLOGICAL

MATERIALS, AND

IMMUNOASSAY

REAGENTS

INNOVATIVE TOOLS
AND TECHNOLOGIES
TO PROVIDE
ASSURANCE IN
DIAGNOSTIC ASSAY
PERFORMANCE AND
TEST RESULTS

FOR MORE
INFORMATION, PLEASE
VISIT OUR WEBSITE:
WWW.SERACARE.COM