

APPLICATION NOTE

Comparison of KPL BacTrace® Anti-Salmonella CSA-Plus Antibody to Two Other Anti-Salmonella species Antibodies in an Indirect ELISA.

Yves Konigshofer, PhD, Alice Ku, Farol L. Tomson, PhD and Seth B. Harkins, PhD

Food poisoning due to *Salmonella* contamination of undercooked meat, usually chicken and eggs, is a major concern for human health and agribusiness. Infection with *Salmonella* causes Salmonellosis characterized by diarrhea, fever and abdominal cramps 12-72 hours after ingestion. According to the US Centers for Disease Control (CDC) the annual incidence of Salmonellosis in the US is approximately 42,000 cases; however, many infections are not reported and the actual rate is thought to be ~1.2 million cases.

There are over 2500+ *Salmonella* serovars and all can cause Salmonellosis. Most (99.5%) of the *Salmonella* serovars are from one subspecies, *Salmonella enterica* subsp. *enterica*. To classify the *Salmonella* serovars further, the serovars are grouped by the O-antigen(s) they express. Historically, O-antigens were designated by letters; however, with the discovery of over 20 O-antigens, the nomenclature has changed to designate O-antigens by number (Grimont). O-antigen expression is important because the O-antigen dictates the host antibody response. Thus, for *Salmonella* detection purposes, researchers need multiple antibodies to detect the various serovars.

In order to eliminate the blending of multiple antibodies specific to every *Salmonella* O-antigen, SeraCare has developed a polyclonal antibody that recognizes a majority of the *Salmonella* O-antigens: BacTrace® Anti-Salmonella CSA-Plus (Cat. No. 01-91-90). The antibody was developed in a two pronged approach. First, goats were immunized using multiple serovars of *Salmonella* that represent the various O-antigens. The broad coverage of O-antigens in the immunogen ensures that no one antigen becomes immuno-dominant. Second, the antibody is purified using a series of affinity purifications utilizing proprietary Encapsulated Column Affinity Purification (ECAP) technology. Utilization of ECAP reduces cross-reactivity while maintaining potency. Combined, these techniques provide Anti-Salmonella CSA-Plus with several advantages over other anti-Salmonella antibodies: 1) Reactivity to most *Salmonella* O-antigens; and 2) High specificity to *Salmonella* with low cross-reactivity to non-*Salmonella* bacteria due to SeraCare's affinity purification methods.

To help researchers visualize the advantages of using BacTrace antibodies in their applications, we have compared Anti-Salmonella CSA-Plus to two commercially available Anti-Salmonella species antibodies in an indirect ELISA format. One of the commercial antibodies is a rabbit polyclonal IgG while

the other is a mouse monoclonal IgG2a with specificity to the following O-antigens: A, B, C, D, E, F & G Groups. Assays were performed to compare both the sensitivity and specificity of the antibodies.

MATERIALS AND METHODS

1. ELISA
 - a. Nunc-Immuno Maxisorp high binding microwell plates (8 well strip format) were coated with 100 µL of various heat-killed bacteria in Phosphate Buffered Saline (PBS, pH 7.4) normalized to OD 7.0 at 650 nm and blocked with 1% Bovine Serum Albumin (BSA) (1:10 dilution of KPL's 10% BSA Diluent/Blocking Solution (Cat# 50-61-00)). Each ELISA plate contained an eight well strip used as a normalization control that was coated with 100 µL of KPL BacTrace *Salmonella* Typhimurium Positive Control (Cat# 50-74-01) and blocked with 1% BSA. Specific bacterial strains used in plate coating are listed in Table 1.
 - b. Primary antibodies: KPL BacTrace (Goat) Anti-Salmonella CSA-Plus (Cat# 01-91-90), USBiological (Rabbit) Anti-Salmonella species (Cat# S0060-29) and Abnova (Mouse) *Salmonella* species (A, B, C, D, E, F & G Groups) monoclonal antibody, clone B343M (Cat# MAB4900). The primary antibodies were diluted to 2 µg/mL in 1% BSA.
 - c. Either 100 µL of primary antibody was added to wells (n=5) or 100 µL of 1% BSA was added to wells (n=3) as a negative background control. Wells were incubated for 1 hour at room temperature.
 - d. Wells were washed using a 1X washing solution (1:20 dilution of KPL's Wash Solution (Cat# 50-63-00)) in a BioTek 405 TS plate washer. The wash cycle includes 3 automatic washes followed by a 9 minute incubation in washing solution ending with 3 more automatic washes.
 - e. 100 µL of the appropriate secondary antibody (KPL's Anti-Goat IgG (H+L) Antibody, Peroxidase Labeled (Cat# 14-13-06), KPL's Anti-Rabbit IgG (H+L) Antibody, Peroxidase Labeled or KPL's Anti-Mouse IgG (H+L) Antibody, Human Serum Adsorbed and

Peroxidase Labeled (Cat# 074-1806) was added to wells in 1% BSA and incubated for 30 minutes.

- f. Wells were washed as in step d.
- g. 100 μ L of KPL's ABTS ELISA HRP substrate (Cat# 50-66-18) was added to each well. Development time varied depending on the primary and secondary antibody combination. Development was stopped with 100 μ L of 1X ABTS Stop Solution (1:5 dilution of KPL's ABTS Peroxidase Stop Solution (Cat# 50-85-01)).
- h. Plates were read using a Molecular Devices Versamax Tunable Plate Reader set at 405 nm.
- i. Individual absorbance readings were background subtracted before averaging. Error bars represent the standard deviation. Unless noted otherwise, each ELISA plate was normalized so that the absorbance value of the *S. Typhimurium* control strip was 2.

RESULTS AND DISCUSSION

Sensitivity Comparison

1. We developed a broad-spectrum anti-*Salmonella* antibody which recognizes the numerous O-antigens expressed by *Salmonella*. The sensitivity of this antibody was measured in an indirect ELISA and compared with two commercially available antibodies. Figure 1 shows the comparison of the normalized absorbance values obtained using various heat-killed *Salmonella* serovars as antigens. The results show that KPL Anti-*Salmonella* CSA-Plus antibody has greater sensitivity to most serovars compared to a rabbit polyclonal. In particular, the KPL Anti-*Salmonella* CSA-Plus showed more sensitivity to O-Groups 8 (C²-C³); 9 (D¹); 3, 10 (E¹); 1, 3, 19 (E⁴) and 13 (G). While the rabbit polyclonal was generally less sensitive than KPL Anti-*Salmonella* CSA-Plus, the mouse monoclonal was very sensitive to a particular strain or not sensitive at all. The mouse monoclonal recognized O-Groups 9 (D¹); 3, 10 (E¹) and 1, 3, 19 (E⁴) well. However, even within the same O-Group, 4 (B) and 7 (C¹), the mouse monoclonal had mixed sensitivity results. Polyclonal antibodies are a better choice for researchers needing broad spectrum bacterial detection because they recognize more epitopes than monoclonal antibodies.

Specificity Comparison

Researchers also value low antibody cross-reactivity. Polyclonal antibodies generally have more cross-reactivity than monoclonal antibodies, and the cross-reactivity can be detrimental when it leads to false positives. However, the cross-reactivity of a polyclonal can be significantly reduced when antigen-specific affinity purification methods are used (i.e., ECAP technology). In contrast to the KPL BacTrace antibody line, most anti-bacterial polyclonal antibodies are purified using Protein G or equivalent. While the resulting antibody is "affinity purified", the polyclonal antibody still contains non-specific antibody which can cause cross-reactivity in an immuno-assay.

To demonstrate the effect of ECAP on reducing cross-reactivity, the KPL Anti-*Salmonella* CSA-Plus antibody was compared to a mouse monoclonal Anti-*Salmonella* species antibody and a rabbit polyclonal Anti-*Salmonella* species antibody, both of which were presumably purified by Protein G. Figure 2 compares the normalized absorbance values obtained in an indirect ELISA with possible cross-reacting bacterial antigens (heat-killed bacterial strains). It is evident in Figure 2 that the rabbit polyclonal shows higher levels of cross-reactivity than either the polyclonal Anti-*Salmonella* CSA-Plus or the mouse monoclonal antibody. Indeed the level of background exhibited by the rabbit polyclonal antibody would obscure many of the positive signals shown in Figure 1. In contrast, the KPL BacTrace Anti-*Salmonella* CSA-Plus shows cross-reactivity which is almost equal to that of a monoclonal antibody.

The results demonstrate that careful antibody development, from immunogen selection through antigen specific affinity purification, produces higher quality antibodies. KPL Anti-*Salmonella* CSA-Plus reacts with higher sensitivity to *Salmonella* species than either a rabbit polyclonal or a mouse monoclonal antibody to a broad range of *Salmonella* serovars. In addition, Anti-*Salmonella* CSA-Plus has the highly desirable trait of low cross-reactivity bordering on the levels typically observed for monoclonal antibodies.

Table 1

Strain/Serovar	ATCC	O-Antigen#
<i>S. Arizonae</i>	13314	51 [IIIa, <i>Salmonella enterica</i> subsp. <i>arizonae</i>]
<i>S. Maartensdijk</i>	15790	40 (R) [IIIa, <i>Salmonella enterica</i> subsp. <i>arizonae</i>]
<i>S. Diarizonae</i>	29934	[IIIb, <i>Salmonella enterica</i> subsp. <i>diarizonae</i>]
<i>S. Harmelen</i>	15783	51 [IV, <i>Salmonella enterica</i> subsp. <i>houtenae</i>]
<i>S. Ochsenzoll</i>	29932	16 (I) [IV, <i>Salmonella enterica</i> subsp. <i>houtenae</i>]
<i>S. Paratyphi A</i>	9150	2 (A)
<i>S. Newington</i>	29628	4 (B)
<i>S. Sloterdijk</i>	15791	4 (B)
<i>S. Typhimurium</i>	14028	4 (B)
<i>S. Choleraesuis</i>	10708	7 (C ₁)
<i>S. Infantis</i>	51741	7 (C ₁)
<i>S. Tennessee</i>	10722	7 (C ₁)
<i>S. Hadar</i>	51956	8 (C ₂ -C ₃)
<i>S. Kentucky</i>	9263	8 (C ₂ -C ₃)
<i>S. Muenchen</i>	8388	8 (C ₂ -C ₃)
<i>S. Tallahassee</i>	12002	8 (C ₂ -C ₃)
<i>S. Newport</i>	6962	8 (C ₂ -C ₃)
<i>S. Berta</i>	8392	9 (D ₁)
<i>S. Enteritidis</i>	49214	9 (D ₁)
<i>S. Gallinarum</i>	9184	9 (D ₁)
<i>S. Typhi</i>	9992v	9 (D ₁)
<i>S. Pullorum</i>	9120	9 (D ₁)
<i>S. Give</i>	9268	3, 10 (E ₁)
<i>S. Newington</i>	29628	3, 10 (E ₁)
<i>S. Senftenberg</i>	8400	1, 3, 19 (E ₄)
<i>S. Simsbury</i>	12004	1, 3, 19 (E ₄)
<i>S. Rubislaw</i>	10717	11 (F)
<i>S. Cubana</i>	12007	13 (G)
<i>S. Putten</i>	15787	13 (G)
<i>S. Havana</i>	NCTC 6086*	13 (G)
<i>S. Mississippi</i>	NCTC 6487*	13 (G)
<i>S. Florida</i>	10727	6, 14 (H)
<i>S. Minnesota</i>	9700	21 (L)
<i>Citrobacter braakii</i>	6750	N/A
<i>Citrobacter freundii</i>	6879	N/A
<i>Enterobacter aerogenes</i>	13048	N/A
<i>Enterobacter cloacae</i>	13047	N/A
<i>Escherichia coli</i>	25922	N/A
<i>E. coli</i> O157	35150	N/A
<i>Hafnia alvei</i>	29926	N/A
<i>Klebsiella pneumoniae</i>	27736	N/A
<i>Kluyvera ascorbata</i>	33433	N/A
<i>Proteus hauseri</i>	13315	N/A
<i>Serratia marcescens</i>	14756	N/A
<i>Shigella boydii</i>	8702	N/A
<i>Shigella flexneri</i>	12025	N/A
<i>Shigella sonnei</i>	11060	N/A
<i>Staphylococcus aureus</i>	6538	N/A
<i>Staphylococcus epidermidis</i>	12228	N/A
<i>Yersinia ruckeri</i>	29473	N/A

*National Collection of Type Cultures, Public Health England #Based on WHO reference, format: WHO numeric O-group (historic letter O-group). Unless otherwise noted *Salmonella* serovars are *Salmonella enterica* subsp. *enterica*.




Figure 1: Normalized absorbance data for an indirect ELISA using KPL Anti-Salmonella CSA-Plus (green), a rabbit polyclonal Anti-Salmonella species (blue) and a mouse monoclonal Anti-Salmonella species (purple) against various Salmonella serovars. Roman numbers represent different Salmonella enterica subspecies. Numbers and letters in parenthesis represent different conventions of Salmonella enterica subsp. enterica O-antigens. Data and error bars represent the mean \pm 1 SD (n=5).




Figure 2: Cross-reactivity data showing normalized absorbance data for an indirect ELISA using KPL Anti-Salmonella CSA-Plus (green), a rabbit polyclonal Anti-Salmonella species (blue) and a mouse monoclonal Anti-Salmonella species (purple) against various bacterial strains. Data and error bars represent the mean \pm 1 SD (n=5).

REFERENCES

Grimont, PAD and Weill, FX. "Antigenic Formulae of the Salmonella Serovars, 9th edition" 2007. WHO Collaborating Centre for Reference and Research on *Salmonella*, World Health Organization, Institut Pasteur, Paris, France.

RELATED PRODUCTS

Product	Cat. No.
Anti-Salmonella CSA+ Plus Antibody	01-91-90
Anti-Salmonella CSA-1 Antibody	01-91-99
Anti-Salmonella CSA-1 Magnetic Beads	082-01-91-99
Anti-Salmonella CSA-1 Latex Beads	082-02-91-99
Salmonella Typhimurium Positive Control	50-74-01



A SeraCare Company

www.kpl.com
800.638.3167
301.948.7755