

AccuSpan™ BKV

Linearity Panel

2410-0344 / Batch #10688243

OVERVIEW


The AccuSpan™ BKV Linearity Panel 2410-0344 / Batch #10688243 is a seven-member panel consisting of six members representing serial log dilutions of cultured BK virus, with established reactivity for BKV DNA, in BKV DNA negative diluent. This panel also consists of one negative member prepared from the diluent. The diluent was prepared from normal human plasma that was 0.2 micron filtered. Sodium azide (0.09%) was added as a preservative.

Roche cobas® 6800/8800 BKV Assay results are reported for each panel member. Linearity is shown graphically by plotting observed results against expected results. The WHO International Standard was tested in the same run as the AccuSpan™ BKV Linearity Panel members. Both expected and observed results for the standards are reported; the expected values from the WHO standard are based upon application of a dilution factor to the WHO assigned value.

For Research Use Only. Not for use in diagnostic procedures. Data are offered for informational purposes. LGC SeraCare Life Sciences does not claim that others can duplicate test results exactly.

CAUTION: These materials have not been treated and should be considered biohazardous. Follow Universal Precautions.

AccuSpan™ BKV Linearity Panel



BKV DNA results were obtained using the Roche cobas® 6800/8800 BKV test method. A line of best fit is shown.

AccuSpan™ BKV

Linearity Panel

2410-0344 / Batch #10688243

BKV DNA

 Roche cobas® 6800/8800
 BKV Test¹

Panel Member	Results (IU/mL)	log (IU/mL)
01	2.63 x 10⁷	7.42
02	1.97 x 10⁶	6.29
03	2.07 x 10⁵	5.32
04	1.64 x 10⁴	4.21
05	1.84 x 10³	3.27
06	2.10 x 10²	2.32
07	TND	TND
Test Date	07-Feb-2024	
Test Site	RL	
Test Kit Range	21.5 – 100,000,000 IU/mL 1.33 – 8.00 log IU/mL	
Kit Part Code	NA	
Kit Lot No.	NA	
Kit Exp. Date	NA	

¹Results are reported as the mean result of three replicates; positive/reactive results are noted in bold red.
 TND = Target Not Detected; RL = Reference Lab; NA = Not Available

AccuSpan™ BKV

Linearity Panel

2410-0344 / Batch #10688243

WHO International Standard: 1ST BKV International Standard (NIBSC Code: 14/212)

Sample ID	Expected Values (log IU/mL)	Observed Values on Roche cobas® 6800/8800	
		BKV Test (log IU/mL) ¹	% Difference ²
01	6.18	6.20	0.37
02	5.18	5.18	0.11
03	4.18	4.19	0.34
Test Date		07-Feb-2024	
Test Site		RL	
Test Kit Range		21.5 – 100,000,000 IU/mL 1.33 – 8.00 log IU/mL	
Kit Part Code		NA	
Kit Lot No.		NA	
Kit Exp. Date		NA	

¹WHO panel was tested in the same test run as the AccuSpan™ BKV Linearity Panel members. Samples were run in triplicate.

Positive/reactive results are noted in bold red.

²Percentage difference is how much the observed log concentration differs from the expected log concentration. Values calculated for reference only. Laboratories may use the data to apply a correction factor to the test results.

RL = Reference Lab; NA = Not Available

The package insert for this panel can be found at www.seracare.com.

A printed copy of the package insert or data sheet may be requested by email at CDx-Info@LGCGroup.com or by phone at 508.244.6400.