



Refractrol-SP Serum Protein Assay Sheet

Product Description: Refractrol-SP Normal
Material Number: 84653
Batch: 10735157
Expiration Date: 21 MAR 2028
Storage Temperature: Ambient Room Temperature (typically 15-30 °C)

Assay values obtained for Reichert Digital and Analog models temperature compensated:

Average Concentration	Uncertainty	Units
7.8	+/- 0.4	g/dL

QC Testing:

Assay values obtained for Atago models when operated at room temperature:

Average Concentration	Uncertainty	Units
7.5	+/- 0.6	g/dL

Safety Precautions: Use Universal Precautions when handling this material.

Approval:

[Signature] 19 JUN 2025
Prepared By Date

[Signature] 19 JUN 2025
QA Verified By Date