

Refractrol-SP Serum Protein Assay Sheet

Product Description: Refractrol-SP High
 Material Number: 84654
 Batch: 10735158
 Expiration Date: 21 MAR 2028
 Storage Temperature: Ambient Room Temperature (typically 15-30 °C)

Assay values obtained for Reichert Digital and Analog models temperature compensated:

Average Concentration	Uncertainty	Units
10.3	+/- 0.4	g/dL

QC Testing:

Assay values obtained for Atago models when operated at room temperature:

Average Concentration	Uncertainty	Units
10.4	+/- 0.6	g/dL

Safety Precautions: Use Universal Precautions when handling this material.

Approval:

[Signature] _____ Date *19 JUN 2025*

[Signature] _____ Date *19 Jun 2025*