

Onco-Ref™ FFPE Cell Slide Reference Standard BRCA1 c.3481_3491del11 p.E1161Fs OMIM: 113705.0025

Catalog Number:	ASE-6112F
Description:	The Onco-Ref™ FFPE Cell Slides is a biologically-relevant molecular reference standard that can be directly incorporated into your sample processing workflows to optimize protocols, evaluate assay sensitivity and specificity, and analyze the impact of workflow changes on downstream analysis. Our FFPE reference standards allow you to optimize your sample processing protocols for formalin-fixed specimens, without the need to use precious clinical samples.
Gene Name:	BRCA1
Mutation Id:	OMIM: 113705.0025
Nucleotide Change:	c.3481_3491del11
Amino Acid Change:	p.E1161Fs
Frequency:	Heterozygous
Size:	1 pack (20 slides)
Format:	FFPE
Volume/Dimension:	5 µm thick section with 7.0mm diameter cell core
Quantity:	n/a
Concentration:	n/a
Cell(s):	HCT116
Shelf Life:	36 months
Shipping Temperature:	Ambient
Fixation Process:	10% buffered formalin solution (4% formaldehyde); 24 hrs; room temperature
Storage Temperature:	-20°C
Quality Control:	Genotype: Sanger sequencing of locus-specific PCR (cell line) Quality: > 60% cell density (H&E staining of 5µm cell slide) Manufacturing: ISO 9001:2015 and ISO 13485:2016 Certified
Regulatory:	Research Use Only. Not Intended for human or animal diagnostic or therapeutic use.
Restrictions:	This is a limited use product. This product, any material that contains this product in whole or in part, any progeny, modification or derivative of this product, any cell or animal made or modified by using this product, and any progeny, modification or derivative of such cell or animal may not be transferred by the purchaser to any other person, entity or any of the above to perform services for the benefit of any other person or entity or for commercial purpose of any kind. This product may only be used by the purchaser for its internal research for use as a research tool for research purposes.

DS-ASE-6112F-V3