

Onco-Ref™ Genomic DNA Reference Standard RKO WT

Catalog Number:	ASO-6051-1
Size:	1 vial
Format:	gDNA
Volume/Dimension:	100 µl
Quantity:	5.0 µg
Concentration:	50 ng/µl
Cell(s):	RKO
Shelf Life:	24 Months
Buffer:	Tris-EDTA (10 mM Tris-HCl, 1 mM EDTA), pH 8.0
Shipping Temperature:	Ambient
Storage Temperature:	-20°C
Quality Control:	Genotype: Whole exon sequencing Quality: Gel Electrophoresis Quantification: Qubit 4.0 fluorometer Manufacturing: ISO 9001:2015 and ISO 13485:2016 Certified
Freeze/Thaw:	No more than 3 cycles
Regulatory:	Research Use Only. Not Intended for human or animal diagnostic or therapeutic use.
Restrictions:	This is a limited use product. This product, any material that contains this product in whole or in part, any progeny, modification or derivative of this product, any cell or animal made or modified by using this product, and any progeny, modification or derivative of such cell or animal may not be transferred by the purchaser to any other person, entity or any of the above to perform services for the benefit of any other person or entity or for commercial purpose of any kind. This product may only be used by the purchaser for its internal research for use as a research tool for research purposes.