KPL ReserveAP™ Phosphatase-Labeled Affinity Purified Antibody To Mouse IgG (H+L) (Human Serum Adsorbed)

Produced in Goat

<u>Catalog No.</u> <u>Size</u> 5220-0468 (4751-1806) <u>Size</u> 1.0 mg/mL



Affinity purified antibody isolated from a pool of serum from goats immunized with purified mouse IgG was labeled with phosphatase using a modified glutaraldehyde procedure (1).

FORM/STORAGE

Liquid. Store refrigerated at 2 - 8°C. Stable for a minimum of 1 year when stored at 2 - 8°C.

STABILIZER AND PRESERVATIVE

Bovine serum albumin (BSA) is added as a protein stabilizer. A proprietary anti-bacterial agent is added to preserve the product. Non-sterile.

ANTIBODY CONCENTRATION

The concentration of affinity purified antibody is 1 mg as determined by UV absorbance at 280 nm.

E/P RATIO

Molar enzyme/antibody protein ratio is from 2.2:1 to 2.9:1.

SPECIFICITY/CROSS REACTIVITY

Tested by gel diffusion and ELISA techniques as applicable. This product reacts specifically with mouse IgG and may recognize other immunoglobulin types that have light chains in common with IgG. Reactivity to IgG subclasses has not been tested. Antibodies to mouse IgG may cross-react with immunoglobulins of other mammalian species if common binding sites are shared. Cross-reactivity with human serum has been minimized with affinity procedures.

SUGGESTED WORKING DILUTIONS

Optimal working concentration should be determined experimentally. Prepare working dilution in TBS or other buffer such as BSA or KPL Milk Diluent/Blocking Solution (See RELATED PRODUCTS) immediately before use. These buffers are not recommended for long term storage. Suggested starting dilutions are as follows.

Note: In many cases, the antibody may be diluted further than indicated.



ELISA:

1/1000 – 1/5000 1.0 μg/mL - 0.2 μg/mL

Blotting:

1/5000- 1/20000 0.2 μg/mL - 0.05 μg/mL

Histo/Cytochemical Procedures:

1/200-1/500 5.0 μg/mL - 2.0 μg/mL

PRODUCT SAFETY AND HANDLING

This product is considered non-hazardous as defined by The Hazard Communication Standard (29 CFR 1910.1200). Avoid contact with skin and eyes. In case of contact or spillage, clean with copious amounts of water. Disposal via sanitary sewer.

REFERENCES

1. Voller, A., et al. (1976). <u>A Manual of Clinical</u> <u>Immunology.</u> American Society for Microbiology, 509-510.

RELATED PRODUCTS	CAT No.
KPL 10% BSA	5140-0006 (50-61-00)
Diluent/Blocking Solution	
Concentrate	
KPL Coating Solution	5150-0014 (50-84-00)
Concentrate	
KPL 20X Wash Solution	5150-0008 (50-63-00)
Concentrate	
KPL SureBlue 1-Component	5120-0076 (52-00-02)
TMB Microwell Substrate	
KPL BluePhos Microwell	5120-0059 (50-88-00)
Substrate	, ,