

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq Breast CNV Mix, +6 copies

Material No: 0710-0412                      Batch Number: 10311353

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Date of Manufacture: 30Jan2018                      Expiration Date: 30Jan2020

Test Method: Qubit dsDNA BR Assay, BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System.

Concentration:      Nominal value: 10 ng/μL  
                             Average measured value: 12.13 ng/μL  
 Fill volume:        20 μL

Gene	Number of Additional Copies by Digital PCR
ERBB2	6.50
FGFR3	6.25
MYC	6.50

Approval:



Prepared By

02/05/2018

Date