

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq Breast CNV Mix, +12 copies

Material No: 0710-0413 Batch Number: 103111354

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Date of Manufacture: 30Jan2018 Expiration Date: 30Jan2020

Test Method: Qubit dsDNA BR Assay, BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System.

Concentration: Nominal value: 10 ng/μL
 Average measured value: 11.83 ng/μL
 Fill volume: 20 μL

Gene	Number of Additional Copies by Digital PCR
ERBB2	12.53
FGFR3	12.15
MYC	12.20

Approval:



Prepared By

02/05/2018

Date