

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® ctDNA EGFR G719S Mutation Mix AF1%

Material No: 0710-0712                      Batch Number: 10396324

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and a synthetic DNA construct

Date of Manufacture: 11-Feb-2019                      Expiration Date: 11-Feb-2021

Concentration (Qubit dsDNA BR Assay): Nominal value: 10 ng/μL; Average measured value: 14.8 ng/μL

Volume: 20 μL

Storage: -20 °C

Digital PCR testing using BioRad QX200™ Droplet Digital™ PCR System:

| Gene ID | COSMIC Identifier | Amino Acid Change | Average AF% |
|---------|-------------------|-------------------|-------------|
| EGFR    | COSM6252          | G719S             | 1.00        |

Approval:



Prepared By

12-Feb-2019

Date