

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 21, 18 and 13 Multi-Analyte Aneuploidy Reference Material

Material No: 0720-0002

Batch Number: 10401901

Material Description: Mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Date of Manufacture: 14-Mar-2019

Expiration Date: 10-Apr-2021

Concentration (Qubit dsDNA
BR Assay): 25.1 ng/mL

Volume: 1.0 mL

Fetal Fraction (In house
digital PCR assays run on
BioRad QX-200 system): 33.5% (T13: 10.1 %, T18: 13.2 %, T21: 15.6 %)

Trisomy Call for
Chromosome 21, 18 and 13
by NGS: Pass

Storage: 2-8 °C

Approval:

A handwritten signature in blue ink, appearing to read "J. Blad".

15-Mar-2019

Prepared By

Date