

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 18 Aneuploidy Reference Material

Material No: 0720-0018

Batch Number: 10395757

Material Description: Mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Date of Manufacture: 25-Jan-2019

Expiration Date: 30-Nov-2022

Concentration (Qubit dsDNA
BR Assay): 29.9 ng/mL

Volume: 1.0 mL

Fetal Fraction (In house
digital PCR assays run on
BioRad QX-200 system): 9 %

Trisomy Call for
Chromosome 18 by NGS: Pass

Storage: 2-8 °C

Approval:

A handwritten signature in black ink, appearing to be "J. K. D.", written over a horizontal line.

Prepared By

29-Jan-2019

Date