

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 21 Aneuploidy Reference Material

Material No: 0720-0019

Batch Number: 10387621

Material Description: Mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Date of Manufacture: 17-Dec-2018

Expiration Date: 03-Dec-2022

Concentration (Qubit dsDNA
BR Assay): 26.7 ng/mL

Volume: 1.0 mL

Fetal Fraction (In house
digital PCR assays run on
BioRad QX-200 system): 9%

Trisomy Call for
Chromosome 21 by NGS: Pass

Storage: 2-8 °C

Approval:

A handwritten signature in black ink, appearing to read "J. Blum", is written over a horizontal line.

Prepared By

16-Jan-2019

Date