

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 21 Aneuploidy Reference Material

Material No: 0720-0019

Batch Number: 10404997

Material Description: Mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Date of Manufacture: 12-Mar-2019

Expiration Date: 17-Jan-2022

Concentration (Qubit dsDNA
BR Assay): 32.7 ng/mL

Volume: 1.0 mL

Fetal Fraction (In house
digital PCR assays run on
BioRad QX-200 system): 10%

Trisomy Call for
Chromosome 21 by NGS: Pass

Storage: 2-8 °C

Approval:

A handwritten signature in blue ink, appearing to read "SeraCare".

Prepared By

15-Mar-2019

Date