

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq[®] Aneuploidy Negative (Euploid) Reference Material

Material No: 0720-0020

Batch Number: 10245834

Material Description: Mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Date of Manufacture: 09-Feb-17

Expiration Date: 09-Feb-21

Concentration (Qubit dsDNA
BR Assay): 27.6 ng/mL

Volume: 1.0 mL

Fetal Fraction (In house
digital PCR assays run on
BioRad QX-200 system): 12%

Euploid Calls for
Chromosome 13, 18 and 21
by NGS: Pass

Storage: 2-8 °C

Approval:

A handwritten signature in black ink, appearing to be "JBL" followed by a long horizontal stroke.

Prepared By

11-Jan-19

Date