

Revised Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Aneuploidy Negative (Euploid) Reference Material

Material Number: 0720-0020 Batch Number: 10388099

Material Description: A mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Concentration
(Qubit dsDNA BR Assay): 26.7 ng/mL

Volume: 1.0 mL


Date of Manufacture: 11-JAN-2019 Expiration Date: 11-JAN-2024

Fetal Fraction
(In house digital PCR
assays run on a BioRad
QX-200 system): 11 %

Euploid Calls for
Chromosome 13, 18
and 21 by NGS: Pass

Storage: 2-8°C

Approval:

A handwritten signature in black ink, appearing to read "M H I".

09 JUN 2023

Prepared By

Date