

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® Euploid Female - Matched Reference Material  
Material Number: 0720-0170 Batch Number: 10683414  
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).  
Concentration (Qubit dsDNA BR Assay): 31.5 ng/mL  
Volume: 1.0 mL  
Date of Manufacture: 27 OCT 2023 Expiration Date: 27 OCT 2027  
Fetal Fraction by NGS (ChromaPanel): 18.0%  
  
Euploid Call by NGS (ChromaPanel): Pass  
Call for fetal gender by NGS (ChromaPanel): XX  
Storage: 2-8 °C

**Approval:**

Prepared By



09 NOV 2023

Date