

# ACCURUN® Anti-SARS-CoV-2 Controls Kit

**ENSURE ANTIBODY TESTING ACCURACY AND PERFORMANCE**

## ABOUT LGC SERACARE

TRUSTED SUPPLIER  
TO THE DIAGNOSTIC  
TESTING INDUSTRY  
FOR OVER 30 YEARS

HIGH-QUALITY  
CONTROL PRODUCTS,  
RAW BIOLOGICAL  
MATERIALS, AND  
IMMUNOASSAY  
REAGENTS

INNOVATIVE TOOLS  
AND TECHNOLOGIES  
TO PROVIDE  
ASSURANCE IN  
DIAGNOSTIC ASSAY  
PERFORMANCE AND  
TEST RESULTS

FOR MORE  
INFORMATION, PLEASE  
VISIT OUR WEBSITE:  
[WWW.SERACARE.COM](http://WWW.SERACARE.COM)

ACCURUN® Anti-SARS-CoV-2 Controls Kit offers positive materials for antibodies to the SARS-CoV-2 virus. Also included are complementary negative materials. The kit is designed to support assay installation and monitoring of day-to-day assay performance, providing a complete quality solution for SARS-CoV-2 antibody testing.

To support the many SARS-CoV-2 serology assay formats available, the kit is currently offered in two configurations, one designed for IgG-specific assays and the other for total antibody tests. Controls are specifically formulated to exhibit low positive reactivity in true patient-like matrices, to test assay performance near critical signal cutoff boundaries.

Product Description	Series	Material Number	Pack Size	Antibody Coverage
ACCURUN Anti-SARS-CoV-2 Controls Kit	1000	2015-0222	Positive 2 x 3 mL Negative 2 x 3 mL	IgG
	2000	2015-0225		Total

*For In Vitro Diagnostic Use. CE-IVD marked.*

## SARS-COV-2 SEROLOGY SOLUTION:

- **Patient-like Matrices:** Formulated from true human-sourced disease material to mimic patient samples encountered in a daily testing environment
- **Relevant Target Ranges:** Specifically designed to be weakly reactive, low positive reference materials to challenge assay performance
- **Ease of Use:** Liquid-ready format, requiring no reconstitution
- **Long Shelf Life:** 2 year stability at 2 – 8°C
- **Open Vial Stability:** Once opened, stable for 30 days

## LEARN MORE

To learn more about ACCURUN Anti-SARS-CoV-2 Controls and our complete line of SARS-CoV-2 Quality Solutions visit: [www.seracare.com/SARSCoV2](http://www.seracare.com/SARSCoV2) or contact customer service at 800.676.1881.

## ORDERING INFORMATION

To place an order, please contact us at +1.508.244.6400 and 800.676.1881 or email [customerservice@seracare.com](mailto:customerservice@seracare.com)