

ACCURUN® SARS-CoV-2 Antigen Reference Material Kit

ASSAY VERIFICATION AND PERFORMANCE MONITORING FOR SARS-COV-2 ANTIGEN TESTS

ABOUT LGC SERACARE

TRUSTED SUPPLIER
TO THE DIAGNOSTIC
TESTING INDUSTRY
FOR OVER 30 YEARS

HIGH-QUALITY
CONTROL PRODUCTS,
RAW BIOLOGICAL
MATERIALS, AND
IMMUNOASSAY
REAGENTS

INNOVATIVE TOOLS
AND TECHNOLOGIES
TO PROVIDE
ASSURANCE IN
DIAGNOSTIC ASSAY
PERFORMANCE AND
TEST RESULTS

FOR MORE
INFORMATION, PLEASE
VISIT OUR WEBSITE:
WWW.SERACARE.COM

ACCURUN® SARS-CoV-2 Antigen Reference Material Kit is formulated for use with test methods that detect the nucleocapsid (NP) protein of SARS-CoV-2 virus. The kit offers both positive and negative materials for SARS-CoV-2 nucleocapsid antigen tests. The kit is designed to support assay installation and monitoring of day-to-day assay performance, providing a complete quality solution for SARS-CoV-2 antigen testing.

Product Description	Material Number	Pack Size	Antigen Coverage
ACCURUN® SARS-CoV-2 Antigen Reference Material Kit	2015-0231	Positive 5 x 3 mL Negative 5 x 3 mL	Nucleocapsid

Not for In Vitro Diagnostic Use. Research Use Only.

SARS-COV-2 ANTIGEN SOLUTION:

- **Nucleocapsid Target:** Formulated with purified, recombinant SARS-CoV-2 nucleocapsid protein in universal transport media
- **Relevant Ranges:** Specifically designed to be weakly reactive, low positive reference material to challenge assay performance
- **Ease of Use:** Liquid-ready format, requiring no reconstitution

LEARN MORE

To learn more about the ACCURUN SARS-CoV-2 Antigen Reference Material Kit and our complete line of SARS-CoV-2 Quality Solutions, visit: www.seracare.com/SARSCoV2 or contact customer service at 800.676.1881.

ORDERING INFORMATION

To place an order, please contact us at +1.508.244.6400 and 800.676.1881 or email customerservice@seracare.com