

ACCURUN® Anti-SARS-CoV-2 Controls Kits

ENSURE ANTIBODY TESTING ACCURACY AND PERFORMANCE

ABOUT LGC SERACARE

TRUSTED SUPPLIER
TO THE DIAGNOSTIC
TESTING INDUSTRY
FOR OVER 30 YEARS

HIGH-QUALITY
CONTROL PRODUCTS,
RAW BIOLOGICAL
MATERIALS, AND
IMMUNOASSAY
REAGENTS

INNOVATIVE TOOLS
AND TECHNOLOGIES
TO PROVIDE
ASSURANCE IN
DIAGNOSTIC ASSAY
PERFORMANCE AND
TEST RESULTS

FOR MORE
INFORMATION, PLEASE
VISIT OUR WEBSITE:
WWW.SERACARE.COM

ACCURUN® Anti-SARS-CoV-2 Controls Kits offer positive materials for antibodies to the SARS-CoV-2 virus. Also included are complementary negative materials. Each kit is designed to support assay installation and monitoring of day-to-day assay performance, providing a complete quality solution for SARS-CoV-2 antibody testing.

To support the many SARS-CoV-2 serology assay formats available, the kit is offered in three configurations for assays that detect IgM, IgG and Total SARS-CoV-2 antibodies, respectively. Controls are specifically formulated to exhibit low positive reactivity in true patient-like matrices, to test assay performance near critical signal cutoff boundaries.

Product Description	Material Number	Pack Size	Antibody Coverage
ACCURUN Anti-SARS-CoV-2 IgM Controls Kit	2015-0232	Positive 2 x 3 mL Negative 2 x 3 mL	IgM
ACCURUN Anti-SARS-CoV-2 Controls Kit, Series 1000	2015-0222		IgG
ACCURUN Anti-SARS-CoV-2 Controls Kit, Series 2000	2015-0225		Total

For In Vitro Diagnostic Use. CE-IVD marked.

SARS-COV-2 SEROLOGY SOLUTION:

- **Patient-like Matrices:** Formulated from true human-sourced disease material to mimic patient samples encountered in a daily testing environment
- **Relevant Target Ranges:** Specifically designed to be weakly reactive, low positive reference materials to challenge assay performance
- **Ease of Use:** Liquid-ready format, requiring no reconstitution
- **Long Shelf Life:** 2 year stability at 2 - 8°C
- **Open Vial Stability:** Once opened, stable for 60 days

LEARN MORE

To learn more about ACCURUN Anti-SARS-CoV-2 Controls and our complete line of SARS-CoV-2 Quality Solutions visit: www.seracare.com/SARSCoV2 or contact customer service at 800.676.1881.

ORDERING INFORMATION

To place an order, please contact us at +1.508.244.6400 and 800.676.1881 or email CDx-CustomerService@lgcgroup.com