

Product Information Sheet

Accukit™ Myeloid DNA IS

Product Information

Accukit™ Name: Accukit™ Myeloid DNA IS
Catalog No: 2023
Storage Conditions: -20°C
Lot Number: 23042807
Expiration Date: APRIL 28, 2025
Concentration: Mixture of each target provided at 300 copies/μl
Volume: 50 μl per vial

Targeted Regions

1. ABL1 chr9:133748283
2. DNMT3A chr2:25457243
3. FLT3 chr13:28592642
4. FLT3 chr13:28608254
5. FLT3 chr13:28608322
6. IDH1 chr2:209113113
7. IDH2 chr15:90631838
8. JAK2 chr9:5073770
9. JAK2 chr9:5070021
10. KRAS chr12:25398281
11. NRAS chr1:115252204
12. NRAS chr1:115256529
13. TET2 chr4:106164914
14. TP53 chr17:7577559
15. NPM1 chr5:170837547



Date: 4/28/2023

QA

FOR LABORATORY USE ONLY