

## Product Information Sheet

### Accukit™ ctDNA MRD IS

#### Product Information

Accukit™ Name: Accukit™ ctDNA MRD IS  
Catalog No: 2091  
Storage Conditions: -20°C  
Lot Number: 23122002  
Expiration Date: DEC 20, 2025  
Concentration: Mixture of each target provided at 300 copies/μl  
Volume: 10 μl per vial

Target concentration confirmed by dPCR  
Fragmented size confirmed to be 170 bp ± 40 bp by Agilent 2100 Bioanalyzer.

#### Targeted Variants

AKT1\_c.49G>A  
ALK\_c.3522C>A  
ALK\_c.3604G>A  
BRAF\_c.1799T>A  
EGFR\_c.2235\_2249del  
EGFR\_c.2369C>T  
EGFR\_c.2573T>G  
KRAS\_c.34G>T  
NRAS\_c.182A>G  
PIK3CA\_c.3140A>G




Figure 1. Agilent Bioanalyzer trace of fragmented DNA in Accukit ctDNA MRD

*Christoph J. Hahn* Date: 3/7/24  
QA