

Product Information Sheet

Accukit™ Myeloid DNA IS (fragmented)

Product Information

Accukit™ Name: Accukit™ Myeloid DNA IS (fragmented)
Catalog No: 2140
Storage Conditions: -20°C
Lot Number: 24010901
Expiration Date: JAN 09, 2026
Concentration: Mixture of each target provided at 300 copies/μl
Volume: 100 μl per vial

***NOTE:** targets provided as mechanically sheared DNA with an average length of ~170 bp to mimic liquid biopsy samples

Targeted Regions

1. ABL1 chr9:133748283
2. DNMT3A chr2:25457243
3. FLT3 chr13:28592642
4. FLT3 chr13:28608254
5. FLT3 chr13:28608322
6. IDH1 chr2:209113113
7. IDH2 chr15:90631838
8. JAK2 chr9:5073770
9. JAK2 chr9:5070021
10. KRAS chr12:25398281
11. NRAS chr1:115252204
12. NRAS chr1:115256529
13. TET2 chr4:106164914
14. TP53 chr17:7577559
15. NPM1 chr5:170837547



Date: JAN 9, 2024

QA

FOR LABORATORY USE ONLY