



Quality Control Certificate

Product Description: Optitrol NAT HIV Screening

Catalog Number: NT01032

Lot: 10810224 (Box)

10810223 (Vial)

Expiry Date: 30 APR 2028

Application: For *in vitro* diagnostic use

Storage: -40 °C to -20 °C

Human Viral Marker: The controls are prepared from processed human plasma or serum positive for HIV RNA.

HIV stock materials were inactivated by gamma radiation.

CE 2797: This product is manufactured and distributed in conformity with the Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

Batch Release: The lot has been tested according to valid directions. It meets all required specification criteria and was released for sale.

Approval:

11 JUN 2026

Quality Assurance

Date