

Representative Reactivities Table



**CLINICAL
DIAGNOSTICS**

Optitrol Orange

Positive Control

Version: 1

Date: Mar 31, 2025

Analyte	Manufacturer	Assay Name	Method	Unit	Mean	SD	Samples	Labs	Representative Reactivity
Anti-HCV IgG	Siemens	ADVIA Centaur HCV	ChLIA	Index	3.02	0.07	10	1	positive
		Atellica HCV	ChLIA	Index	2.95	0.07	10	1	positive
Anti-HIV 1 IgG	Siemens	ADVIA Centaur HIV Ag/Ab Combo	ChLIA	Index	2.93	0.04	10	1	positive
		Atellica HIV Ag/Ab Combo	ChLIA	Index	3.33	0.03	10	1	positive
HBsAg	Siemens	ADVIA Centaur HBsAg II	ChLIA	Index	3.77	0.13	10	1	positive
		Atellica HBsAg II	ChLIA	Index	4.15	0.06	10	1	positive
Anti-HBc IgG	Siemens	ADVIA Centaur Anti-HBc	ChLIA	Index	2.13	0.13	10	1	positive
		Atellica Anti-HBc Total	ChLIA	Index	1.81	0.03	10	1	positive
Anti-Treponema pallidum IgG	Siemens	ADVIA Centaur Syphilis	ChLIA	Index	4.88	0.08	10	1	positive

The stated mean values were obtained by testing one product batch. They show possible representative values but no target values. Test results may vary among used instruments, assays, reagent lots, and different laboratory protocols.