

Technical Product Report

Product: AccuPlex™ SARS-CoV-2 Verification Panel

Part Number 0505-0129

Lot Number: 10481707

Panel details:

Member	Nominal Concentration Copies/ml	Number of vials per panel
SARS-CoV-2 Member 1	1.00E+05	1
SARS-CoV-2 Member 2	1.00E+04	1
SARS-CoV-2 Member 3	1.00E+03	1
SARS-CoV-2 Negative	NA	1

Date of Expiration: 24-Mar-2022

Storage: 2-8 °C

Matrix: Viral Transport Buffer

Volume: 3 mL per vial

Safety Precautions: Preparation contains heat treated recombinant sindbis virus. Use the Centers for Disease Control recommended universal precautions for handling the material.

Approval:


 Approved By

3/26/2020
 Date