

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® CNV Breast Mix, +3 copies

Material Number: 0710-0411 Batch Number: 10446808

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Concentration (Qubit dsDNA BR Assay): Nominal value: 10 ng/μL
Average measured value: 13.0 ng/μL

Fill Volume: 20 μL

Date of Manufacture: 10 OCT 2019 Expiration Date: 10 OCT 2022

Test Method: BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System

Gene	Number of Additional Copies
ERBB2	3.68
FGFR3	3.65
MYC	3.77

Approval:



Prepared By

Date

28 JUN 2021