

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® CNV Breast Mix, +3 copies

Material Number: 0710-0411 Batch Number: 10640753

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Concentration (Qubit dsDNA BR Assay): Nominal value: 10 ng/μL
Average measured value: 14.2 ng/μL


Fill Volume: 20 μL

Date of Manufacture: 21 SEP 2022 Expiration Date: 21 SEP 2024

Test Method: BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System

Gene	Number of Additional Copies
ERBB2	2.8
FGFR3	3.8
MYC	6.2

Approval:



Prepared By

26 SEP 2022

Date