

# Technical Product Report

*For Research Use Only; Not for use in Diagnostic Procedures*

Product Description: Seraseq® CNV Breast Mix, +12 copies

Material Number: 0710-0413      Batch Number: 10623365

Material Description: Mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs

Concentration (Qubit dsDNA BR Assay): Nominal value: 10 ng/μL  
Average measured value: 15.5 ng/μL


Fill Volume: 20 μL

Date of Manufacture: 25 MAY 2022      Expiration Date: 25 MAY 2024

Test Method: BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System

Gene	Number of Additional Copies
ERBB2	12.6
FGFR3	11.2
MYC	12.9

Approval:



01-JUN-2022

Prepared By

Date