

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description:	SeraSeq® ctDNA Complete Mutation Mix AF5%		
Material Number:	0710-0528	Batch Number:	10775745
Material Description:	A ctDNA-like mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs		
Fill Volume:	25 µL		
Date of Manufacture:	14 AUG 2025	Expiration Date:	14 AUG 2029
Storage:	-20°C		
Concentration (Qubit dsDNA BR Assay):	Nominal value: 10 ng/µL Average measured value: 15.2 ng/µL		
Average fragment size (Agilent Bioanalyzer DNA 1000 Analysis):	157 bp		
Acceptance criteria for average fragment size:	140 – 200 bp		

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Mutation Mix AF5%

Digital PCR testing using
BioRad QX200™ Droplet
Digital™ PCR System:

Gene ID	COSMIC Identifier	Amino Acid Change	Average AF%
AKT1	COSM33765	p.E17K	4.46
BRAF	COSM476	p.V600E	4.73
EGFR	COSM6224	p.L858R	4.88
EGFR	COSM6240	p.T790M	5.02
ERBB2	COSM20959	p.A775_G776insYVMA	4.19
KIT	COSM1314	p.D816V	4.75
KRAS	COSM521	p.G12D	4.98
NCOA4/RET	NA	Translocation	4.78
NRAS	COSM584	p.Q61R	4.89
PIK3CA	COSM775	p.H1047R	4.87
PIK3CA	COSM12464 ¹	p.N1068fs*4	4.87
EML4-ALK	NA	Translocation	4.67
ALK	COSM144250	p.G1202R	4.91
ALK	COSM28055	p.F1174L	4.91
BRCA1	COSM1383519	p.K654fs*47	4.75
BRCA2	COSM1738242	p.R2645fs*3	4.89
EGFR	COSM12370	p.L747_P753>S	5.47
EGFR	COSM6256	p.S752_I759delSPKANKEI	5.39
EGFR	COSM6223	p.E746_A750delELREA	5.59
KRAS	COSM516	p.G12C	4.77
CD74/ROS1	NA	Translocation	5.04
KRAS	COSM554	p.Q61H	4.53

Gene ID	Average CNV in ctDNA ²	Average Additional Copies (per cell) in ctDNA
ERBB2	7.50	5.50
MET	5.75	3.75
MYC	6.64	4.64

NA = not applicable

¹As of June 2019, this mutation is no longer listed in the COSMIC database.

²Compare to a normal CN of 2.00.

Next Generation Sequencing testing using Archer® Reveal ctDNA™ 28 Kit run on an Illumina® MiSeq™ using v2 (2x150 bp) PE chemistry reagents¹:

Gene ID	COSMIC Identifier	Amino Acid Change	AF%
AKT1	COSM33765	p.E17K	3.69
BRAF	COSM476	p.V600E	4.30
EGFR	COSM6224	p.L858R	4.23
EGFR	COSM6240	p.T790M	4.83
ERBB2	COSM20959	p.A775_G776insYVMA	3.79
KIT	COSM1314	p.D816V	5.15
KRAS	COSM521	p.G12D	4.24
NCOA4/RET	NA	Translocation	NA
NRAS	COSM584	p.Q61R	4.68
PIK3CA	COSM775	p.H1047R	4.65
PIK3CA	COSM12464 ²	p.N1068fs*4	3.88
EML4-ALK	NA	Translocation	NA
ALK	COSM144250	p.G1202R	4.65
ALK	COSM28055	p.F1174L	4.74
BRCA1	COSM1383519	p.K654fs*47	NA
BRCA2	COSM1738242	p.R2645fs*3	NA
EGFR	COSM12370	p.L747_P753>S	4.69
EGFR	COSM6256	p.S752_I759delSPKANKEI	5.06
EGFR	COSM6223	p.E746_A750delELREA	4.83
KRAS	COSM516	p.G12C	4.07
CD74/ROS1	NA	Translocation	NA
KRAS	COSM554	p.Q61H	4.22

Gene ID	CNV in ctDNA ³	Additional Copies (per cell) in ctDNA
ERBB2	6.88	4.88
MET	9.62	7.62
MYC	NA	NA

NA = not applicable; AF% and CNV marked NA were not targeted by the panel.

¹NGS was performed as an orthogonal verification step. Parameters used:

DNA input = 50 ng

of samples / flow cell = 1

of total reads / sample = 90M

Average read depth = 13871X

On-target reads = 86.7%

Q30 score = 89.3%

Analysis = Archer Analysis Suite v6.2.7 (default settings except for: Variant Downstream ROI Size of 150 and Read depth Normalization of 10,000,000)

²As of June 2019, this mutation is no longer listed in the COSMIC database.

³Compare to a normal CN of 2.00.

Note: The MET gene is amplified using two synthetic constructs with a small region of overlap between the constructs (see package insert for genomic coordinates). Assays which target this region of overlap may report higher amplification levels.

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Mutation Mix AF5%

Approval:

Prepared By	<u>Shuang Yang</u>	Date	<u>25 AUG 2025</u>
-------------	--------------------	------	--------------------