

Technical Product Report – REVISED

For Research Use Only; Not for use in Diagnostic Procedures

Product Description:	SeraSeq® ctDNA Complete Reference Material AF 5%		
Material Number:	0710-0669	Batch Number:	10659790
Material Description:	A ctDNA-like mixture of human genomic DNA from the reference cell line, GM24385, and synthetic DNA constructs		
Concentration (Qubit dsDNA BR Assay):	Nominal value: 25 ng/mL; Average measured value after extraction using Qiagen QIAamp Circulating Nucleic Acid Kit: 32.2 ng/mL		
Fill Volume:	5.0 mL		
Date of Manufacture:	05 APR 2023	Expiration Date:	05 APR 2027
Storage:	2-8°C		

Technical Product Report – REVISED

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Reference Material AF 5%

Digital PCR testing using
BioRad QX200™ Droplet
Digital™ PCR System:

Gene ID	COSMIC Identifier	Amino Acid Change	Average AF%
AKT1	COSM33765	p.E17K	5.22
BRAF	COSM476	p.V600E	5.26
EGFR	COSM6224	p.L858R	5.30
EGFR	COSM6240	p.T790M	5.03
ERBB2	COSM20959	p.A775_G776insYVMA	5.12
KIT	COSM1314	p.D816V	5.18
KRAS	COSM521	p.G12D	5.35
NCOA4/RET	NA	Translocation	5.07
NRAS	COSM584	p.Q61R	5.47
PIK3CA	COSM775	p.H1047R	4.98
PIK3CA	COSM12464 ¹	p.N1068fs*4	4.98
EML4-ALK	NA	Translocation	4.58
ALK	COSM144250	p.G1202R	4.88
ALK	COSM28055	p.F1174L	4.88
BRCA1	COSM1383519	p.K654fs*47	4.79
BRCA2	COSM1738242	p.R2645fs*3	4.61
EGFR	COSM12370	p.L747_P753>S	6.02
EGFR	COSM6256	p.S752_I759delSPKANKEI	5.14
EGFR	COSM6223	p.E746_A750delELREA	5.72
KRAS	COSM516	p.G12C	5.28
CD74/ROS1	NA	Translocation	4.88
KRAS	COSM554	p.Q61H	4.89

Gene ID	Average CNV in ctDNA ²	Average Additional Copies (per cell) in ctDNA
ERBB2	7.67	5.67
MET	6.11	4.11
MYC	6.27	4.27

NA = not applicable

¹As of June 2019, this mutation is no longer listed in the COSMIC database.

²Compare to a normal CNV of 2.00.

Technical Product Report – REVISED

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Reference Material AF 5%

Next Generation Sequencing testing using Archer® Reveal ctDNA™ 28 Kit run on an Illumina® MiSeq™ using v2 (2x150 bp) PE chemistry reagents^{1,2}:

Gene ID	COSMIC Identifier	Amino Acid Change	AF%
AKT1	COSM33765	p.E17K	5.00
BRAF	COSM476	p.V600E	5.00
EGFR	COSM6224	p.L858R	4.00
EGFR	COSM6240	p.T790M	4.80
ERBB2	COSM20959	p.A775_G776insYVMA	4.10
KIT	COSM1314	p.D816V	5.90
KRAS	COSM521	p.G12D	5.90
NCOA4/RET	NA	Translocation	NA
NRAS	COSM584	p.Q61R	5.70
PIK3CA	COSM775	p.H1047R	5.30
PIK3CA	COSM12464 ³	p.N1068fs*4	5.20
EML4-ALK	NA	Translocation	NA
ALK	COSM144250	p.G1202R	4.70
ALK	COSM28055	p.F1174L	4.70
BRCA1	COSM1383519	p.K654fs*47	NA
BRCA2	COSM1738242	p.R2645fs*3	NA
EGFR	COSM12370	p.L747_P753>S	5.20
EGFR	COSM6256	p.S752_I759delSPKANKEI	4.50
EGFR	COSM6223	p.E746_A750delELREA	4.50
KRAS	COSM516	p.G12C	5.00
CD74/ROS1	NA	Translocation	NA
KRAS	COSM554	p.Q61H	4.30

Gene ID	CNV in ctDNA ⁴	Additional Copies (per cell) in ctDNA
ERBB2	6.22	4.22
MET	7.12	5.12
MYC	NA	NA

NA = not applicable; AF% and CNV marked NA were not targeted by the panel.

¹NGS was performed as an orthogonal verification step. Parameters used:

DNA input = 50 ng

of samples / flow cell = 7

of total reads / sample = 5.5 M

Average read depth = 10368X

On-target reads = 96.7%

Q30 score = 91.49%

Analysis = Archer Analysis Suite v6.2.7 (default settings except for: N/A)

²Please see the poster from NIST for more information about assay sensitivity:

<https://digital.seracare.com/multilab-assessment-reference-materials-ctdna-poster2018>

³As of June 2019, this mutation is no longer listed in the COSMIC database.

⁴Compare to a normal CNV of 2.00.

Technical Product Report – REVISED

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® ctDNA Complete Reference Material AF 5%

Note: The MET gene is amplified using two synthetic constructs with a small region of overlap between the constructs (see package insert for genomic coordinates). Assays which target this region of overlap may report higher amplification levels.

Approval:

BL

03/29/2024

Prepared By

Date