

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Solid Tumor CNV Mix, +3 copies

Material Number: 0710-2866 Batch Number: 10792795

Material Description: Mixture of human genomic DNA from reference cell line, GM24385, and biosynthetic DNA constructs

Concentration: Nominal value: 10 ng/μL
(Qubit dsDNA BR Assay): Average measured value: 14.1 ng/μL

Fill Volume: 20 μL

Date of Manufacture: 07 NOV 2025 Expiration Date: 07 NOV 2029

Test Method: BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System

Gene	Additional Copies (dPCR)	Additional Copies (NGS)
AKT2	3.27	2.41
BRAF	2.46	2.62
EGFR	3.56	4.07
ERBB2	2.94	3.65
FGFR3	2.84	3.92
KIT	3.03	3.12
KRAS	3.32	2.93
MET	3.83	7.87
MYC	2.80	3.40
MYCN	3.64	3.68
NTRK1	2.98	N/A
PIK3CA	2.61	2.43

Approval:


Prepared By

04 DEC 2025
Date