

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Solid Tumor CNV Mix, +3 copies

Material Number: 0710-2866 Batch Number: 10812880

Material Description: Mixture of human genomic DNA from reference cell line, GM24385, and biosynthetic DNA constructs

Concentration: Nominal value: 10 ng/μL
(Qubit dsDNA BR Assay): Average measured value: 12.4 ng/μL

Fill Volume: 20 μL

Date of Manufacture: 08 MAY 2026 Expiration Date: 08 MAY 2030

Test Method: BioRad PrimePCR Copy Number Assays run on the BioRad QX-200 Droplet Digital PCR System

Gene	Additional Copies (dPCR)	Additional Copies (NGS)
AKT2	2.71	N/A
BRAF	2.22	N/A
EGFR	2.61	3.9
ERBB2	2.57	3.3
FGFR3	2.64	3.2
KIT	2.47	1.0
KRAS	2.37	1.1
MET	2.77	1.9
MYC	2.16	3.5
MYCN	3.19	4.7
NTRK1	2.85	N/A
PIK3CA	2.06	1.2

Approval:



02 JUN 2026

Prepared By

Date