

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 13 Aneuploidy Reference Material

Material Number: 0720-0017 Batch Number: 10510394

Material Description: A mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Concentration
(Qubit dsDNA BR Assay): 25.1 ng/mL

Volume: 1.0 mL

Date of Manufacture: 27 JUL 2020 Expiration Date: 05 JUN 2024

Fetal Fraction
(In house digital PCR
assays run on a BioRad
QX-200 system): 13%

Trisomy Call for
Chromosome 13 by NGS: Pass

Storage: 2-8°C

Approval:

A handwritten signature in blue ink, appearing to be "J. M. ...".

Prepared By

17 AUG 2020
Date