

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 13 Aneuploidy Reference Material

Material Number: 0720-0017 Batch Number: 10662441

Material Description: A mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Concentration
(Qubit dsDNA BR Assay): 25.8 ng/mL

Volume: 1.0 mL

Date of Manufacture: 20 APR 2023 Expiration Date: 20 MAR 2027

Fetal Fraction
(NIPTFY® Focus Plus) 9.8 %

Trisomy Call for
Chromosome 13 by
NIPTFY® Focus Plus: Pass

Storage: 2-8°C

Approval:

A handwritten signature in black ink, appearing to be "M H".

Prepared By

26 APR 2023

Date