

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 18 Aneuploidy Reference Material

Material Number: 0720-0018 Batch Number: 10685006

Material Description: A mixture of human genomic DNA (maternal and fetal) in a commutable matrix (simulated plasma).

Concentration
(Qubit dsDNA BR Assay): 22.4 ng/mL

Volume: 1.0 mL

Date of Manufacture: 10-AUG-2023 Expiration Date: 10-AUG-2027

Fetal Fraction
(In house digital PCR
assays run on a BioRad
QX-200 system): 11.0 %

Trisomy Call for
Chromosome 18 by NGS: Pass

Storage: 2-8°C

Approval:

A handwritten signature in black ink, appearing to be "M H".

Prepared By

17 NOV 2023

Date