

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 21 Female - Matched Reference Material
Material Number: 0720-0168 Batch Number: 10670287
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 31.0 ng/mL
Volume: 1.0 mL
Date of Manufacture: 09 JUN 2023 Expiration Date: 03 FEB 2027
Fetal Fraction by NGS (Illumina VeriSeq): 7.9%
Trisomy Call for Chromosome 21 by NGS (Illumina VeriSeq): Pass
Call for fetal gender by NGS (Illumina VeriSeq): XX
Storage: 2-8°C

Approval:

A handwritten signature in blue ink, appearing to be "M. H.", written over a horizontal line.

19 JUN 2023

Prepared By

Date