


Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 21 Female - Matched Reference Material
Material Number: 0720-0168 Batch Number: 10690466
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 26.7 ng/mL
Volume: 1.0 mL
Date of Manufacture: 21 JAN 2024 Expiration Date: 08 FEB 2027
Fetal Fraction by NGS (Sage™ Prenatal Screen): 9.5%
Trisomy Call for Chromosome 21 by NGS (Sage™ Prenatal Screen): Pass
Call for fetal gender by NGS (Sage™ Prenatal Screen): XX
Storage: 2-8°C

Approval:



19 FEB 2024

Prepared By

Date