

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description:	SeraSeq® Euploid Female - Matched Reference Material		
Material Number:	0720-0170	Batch Number:	10820112
Material Description:	A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).		
Concentration (Qubit dsDNA BR Assay):	29.3 ng/mL		
Volume:	1.0 mL		
Date of Manufacture:	07 AUG 2026	Expiration Date:	07 AUG 2030
Fetal Fraction by NGS (Chroma Panel):	22%		
Euploid Call by NGS (Chroma Panel):	Pass		
Call for fetal gender by NGS (Chroma Panel):	XX		
Storage:	2-8 °C		

Approval:

Prepared By		18 AUG 2026
		Date