

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 18 Female - Matched Reference Material
 Material Number: 0720-0172 Batch Number: 10766381
 Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
 Concentration (Qubit dsDNA BR Assay): 34.4 ng/mL
 Volume: 1.0 mL
 Date of Manufacture: 23 JUN 2025 Expiration Date: 11 SEP 2027
 Fetal Fraction by NGS (Chroma): 9.5%

Trisomy Call for Chromosome 18 by NGS (Chroma): Pass
 Call for fetal gender by NGS (Chroma): XX
 Storage: 2-8°C

Approval:



Prepared By

29 Jul 2025

Date