

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 18 Female - Matched Reference Material
Material Number: 0720-0172 Batch Number: 10806799
Material Description: A ccfDNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
Concentration (Qubit dsDNA BR Assay): 27.6 ng/mL
Volume: 1.0 mL
Date of Manufacture: 24 FEB 2026 Expiration Date: 11 SEP 2027
Fetal Fraction by NGS (Illumina VeriSeq): 9%

Trisomy Call for Chromosome 18 by NGS (Illumina VeriSeq): Pass
Call for fetal gender by NGS (Illumina VeriSeq): XX
Storage: 2-8°C

Approval:

A handwritten signature in blue ink, appearing to be "MM".

19 MAR 2026

Prepared By

Date