

Technical Product Report

For Research Use Only; Not for use in Diagnostic Procedures

Product Description: Seraseq® Trisomy 18 Female - Matched Reference Material
 Material Number: 0720-0172 Batch Number: 10820913
 Material Description: A ccfdNA-like mixture of human genomic DNA (matched maternal and fetal) in a commutable matrix (simulated plasma).
 Concentration (Qubit dsDNA BR Assay): 29.8 ng/mL
 Volume: 1.0 mL
 Date of Manufacture: 24 AUG 2026 Expiration Date: 14 JUL 2030
 Fetal Fraction by NGS (CHROMA Panel): 11%

Trisomy Call for Chromosome 18 by NGS (CHROMA Panel): Pass
 Call for fetal gender by NGS (CHROMA Panel): XX
 Storage: 2-8°C

Approval:

Prepared By



15 SEP 2026

Date